WARNING: Breast implants have been associated with the development of cancer of the immune system called Breast Implant Associated Anaplastic Large Cell Lymphoma (BIA-ALCL). More detailed information can be found in the BIA-ALCL section.

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- Breast augmentation (primary and revision): increases breast size and revision surgery to correct or improve the result of a previous breast augmentation surgery.
- Breast Reconstruction (primary and revision): to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast anomaly, as well as revision surgery to correct or improve the results of a previous breast reconstruction surgery.

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- To enhance the aesthetic appearance in patients dissatisfied with the shape/size of their breasts.
- To correct congenital or acquired breast deformities or breast asymmetry.
- To correct or improve the result of a previous breast augmentation/reconstruction surgery.

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Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® are invasive, long-term implantable medical devices for breast augmentation/reconstruction. All implants are composed of a shell, a patch, and highly viscous and highly elastic Establishment Labs silicone gel formula, ProgressiveGel® PLUS or ProgressiveGel® ULTIMA®. The shell is constructed of successive cross-linked layers of silicone elastomer and a low diffusion barrier layer technology that provides these implants their elasticity and integrity. Motiva Implant Matrix® Sterile Silicone Breast Implants are available with and without microtransponder. The microtransponder is placed inside the Motiva Implant Matrix® Sterile Silicone Breast Implants filler material.

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
Compound	Concentration (µg/unit)
VOC (HS-GC/MS)	
Trimethyl silanol	415
Benzene	1.03
SVOC (GC/MS)	
Benzoic acid	5,81
Caprolactam	53.6
4-Chlorobenzoic acid	189
4-Chlorobenzoic acid, trimethylsilyl ester	32.7
2,4-Dichlorobenzoic acid	328.9
Decamethyl cyclopentasiloxane (D5)	120.3
Dodecamethyl cyclohexasiloxane (D6)	748.1
Tetradecamethyl cycloheptasiloxane (D7)	513.8
Hexadecamethyl cyclooctasiloxane (D8)	165.5
Octadecamethyl cyclononasiloxane (D9)	391.5
2,2,4,4,6,8-Hexamethyl-6,8- diphenyl- cyclotetrasiloxane	7,506
Eicosamethyl cyclodecasiloxane (D10)	1,053
2,2,4,4,6,6,8,10-Octamethyl-8, 10-diphenyl cyclopentasiloxane	19,485
Cyclic polydimethylsiloxane oligomer (combined values)	113,059
Dimethylsiloxanemethylphenylsiloxane copolymer (combined values)	54,717
Siloxane** (combined values)	11,729.8
1,3,5,7-Tetramethyltetraphenyl cyclotetrasiloxane	19,884
2,4,6,8,10-Pentamethyl-2,4,6,8, 10-pentaphenyl-1,3,5,7,9,2,4,6,8, 10-pentaoxapentasiloxane isomer (combined values)	50,794
Polymethylphenylsiloxane oligomer (combined values)	23,974
Octamethyl cyclotetrasiloxane (D4)	79.2
2-Ethyl-hexanol	79.3
Tetracosamethyl cyclodecasiloxane (D12)	804
Linear polydimethylsiloxane oligomer (combined values)	377
NVOC (LC/UV)	
Palmitic acid	158.5
Stearic acid	168.2
Erucamide	43.79
Irganox 245	23.7
NVOC (LC/UV-Vis)	
Siloxane** (combined values)	250,375
Di(2-ethylhexyl)phthalate	9,439
Unknown*** (n = 9)	576.87-1,308
GC/MS, gas chromatography/mass spectrometry; LC/UV, liquid chromatography/ultraviolet; LC/UV-Vis, liquid chromatography/ultraviolet visible; NVOC, non-volatile organic compounds; SVOC, semi-volatile organic compounds; VOC, volatile organic compounds.	

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- Increase the breast size and/or
- Reconstruct the breast to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality or
- Revision procedure to correct or improve the results of previous breast augmentation or reconstruction surgery.

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- Do not allow sharp objects, such as scalpels or needles, to contact the device during the implantation or other surgical procedures.
- Do not immerse the implant in iodine solution. If iodine solution is used in the pocket, make sure that it is rinsed thoroughly with deionized water so that no residual solution remains in the pocket.
- Do not allow the implant to contact cauterization devices.
- Do not alter the implant or attempt to repair or insert a damaged implant.
- Ensure no excessive force is applied to a very small shell area during the insertion of the device throughout the insertion. Instead, apply force over as large an area of the implant as possible during insertion. Excessive force can lead to implant damage by either gel fracture or implant rupture.
- Do not use the periumbilical approach to place the implant.
- Do not place more than one implant per breast pocket.
- Do not treat capsular contracture by closed capsulotomy or forceful external compression, which will likely produce implant damage, rupture, folds, and/or hematoma.
- Do not re-use or re-sterilize any product that has been previously implanted. Breast implants are intended for single use only; therefore, ignoring this indication can cause infection, seroma, capsular contracture, reoperation, unsatisfactory aesthetic results, among other related adverse effects that compromise the patient's health.
- Do not use microwave diathermy in patients with breast implants, as it has been associated with tissue necrosis, skin erosion, and implant extrusion.

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Surgical technique and implant selection: Several surgical techniques can be used to perform a silicone gel-filled breast implant augmentation/reconstruction procedure. Therefore, the surgeon is advised to use his/her clinical judgment to choose the best for the patient, consistent with this product insert.

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Patient Counseling Information

This document and the information for the patient must be carefully reviewed before counseling a patient about Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® and breast augmentation surgery. Doctors must thoroughly read and understand the contents of these documents and make sure that any questions or concerns have been resolved before proceeding with the device's use. Breast implant surgery is an elective procedure, and the patient must understand its potential risks and benefits to make an informed decision. For this reason, the patient should be instructed to read the document called "Breast Augmentation and Reconstruction with Motiva Implants® Information for the Patient", which is available on the website <https://ifu.motiva.health/>. The doctor must discuss with the patients the sections on warnings, contraindications, precautions, essential factors to consider, complications², and all other aspects of the document. The physician should inform the patient about the potential complications and that medical management of serious complications may include additional surgery and explantation.

INFORMED CONSENT

It is surgeon's responsibility to document a successful informed decision-making process, providing his/her "Informed Consent Document", which must be signed by the surgeon, the patient and a witness; this document will become part of the patient's medical record. All possible side effects and complications related to surgery should be advised to patients when considering breast augmentation/reconstruction with silicone gel breast implants.

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The SSCP document could be consulted for additional and updated information on the

Motiva Implant Matrix® on the company website <https://ifu.motiva.health/> and EUDAMED website <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI.

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All patients seeking an elective procedure such as breast augmentation must have realistic expectations that focus on improvement rather than perfection. Ask the patient to openly discuss, before surgery, any history of depression or other mental health disorders.

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Mental Health and Elective Surgery – All patients seeking an elective procedure such as breast augmentation must have realistic expectations that focus on improvement rather than perfection. Ask the patient to openly discuss, before surgery, any history of depression or other mental health disorders.

Surgical Setting and Anesthesia – The Augmentation/Reconstruction surgery is performed in a specialized operating room in a clinic/hospital. General anesthesia is commonly used, and local anesthesia with sedation is also an option. You should be sure to check with your surgeon and with the facility where the surgery will take place to become aware of the tests, pre-surgical examinations, and length of time you need to be without food or your routine medications before the surgical procedure.

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Pain – Most women undergoing augmentation or reconstruction with a mammary implant will experience post-operative breast and/or chest pain. While this pain usually recedes in most women as they heal after surgery, it can become a chronic problem in other women. Hematoma, migration, infection, implants that are too large or capsular contracture can cause chronic pain. Sudden, severe pain may be associated with implant rupture. The surgeon should instruct the patient to immediately report if there is significant pain or if pain persists.

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Lactation Difficulties – Some women who undergo breast augmentation can successfully breastfeed and some cannot. Women who undergo mastectomies and then have breast implant-reconstruction surgeries may not be able to breastfeed on the affected side due to loss of breast tissue and the glands that produce milk.

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The clinical symptoms resulting from a bottomed-out implant include asymmetry, upward-pointing nipples, sagging, palpability, etc. Appropriate surgical planning can mitigate the possible causes of bottoming out. Recommendations include a careful and individual assessment of mammary tissues, careful implant selection, application of risk-minimizing surgical techniques, and adequate breast support

1 Loureno FR, Kikuchi IS, Andreoli Pinto T de J. Silicone Gel Bleed on Breast Implants. Open Biomater J. 2011;3(October):14-17. doi: 10.2174/1876502501103010014

2 Berman B, Maderal A, Raphael B. Keloids and Hypertrophic Scars: Pathophysiology, Classification, and Treatment. Dermatol Surg. 2017;Jan;43 Suppl 1:S3-S18. doi: 10.1097/DSS.0000000000000819. PMID: 27347634.

3 Mañero, Ivan M.D.; Montull, Patricia M.D.; Guisantes, Eva M.D. Bottoming Out: A Simple Technique for Correcting Breast Implant Displacement. Plastic and Reconstructive Surgery. December 2009 - Volume 124 - Issue 6 - p 452e-453e

after surgery. Treatments may vary depending on the severity of the complication and range from a simple sub-mammary fixation to the use of additional supporting materials.

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Capsular Contracture – A capsular contracture pertains to hypertrophic scar tissue investing in a foreign body or surgically implanted device, compromising the aesthetic outcome, resulting in pain, breast deformity, and often necessitating further operations⁵. Detection of breast cancer by mammography may also be challenging. Capsular contracture may be more common following infection, hematoma and seroma, and the chance of it happening may increase over time. Capsular contracture occurs more commonly in patients undergoing revision surgery than in patients undergoing primary implantation surgery. Capsular contracture is the most common complication following implant-based breast surgery and is one of the most common reasons for reoperation.

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Patients should also be advised that additional surgery might be needed in cases where pain and/or firmness are severe (i.e., Baker Grades III or IV) and that capsular contracture may happen again after additional surgeries.

A closed capsulotomy (i.e., external manipulation of the capsule to "pop" the tissue capsule and open it up) is used as a standard procedure for treating capsular contracture; however, most manufacturers, including Establishment Labs, contraindicate it because it can cause implant rupture.

Infection – Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation⁶. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are more challenging to treat than infections in tissue without an implant present. If an infection does not respond to antibiotics, the implant may have to be removed, with replacement occurring only after the infection is resolved.

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Often, seromas are reabsorbed by the body over several weeks, but needle drainage is sometimes needed to remove the fluid⁸. While seromas do not increase breast cancer risk, they sometimes heal with scar tissue or calcifications that can raise a concern about mammograms in the future. Seroma symptoms most often appear a week to 10 days after surgery; the area may feel tender and swollen, with a discrete lump and redness arising within a day or two. Early seroma formation is defined as periprosthetic fluid

accumulation within the first postoperative year, whereas the late form is any moment beyond that time⁹.

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Hematoma – A hematoma is a collection of blood within the breast tissue. Hematomas are one of the several complications that may follow a breast-augmentation procedure. Symptoms of hematomas generally include swelling, bruising and pain around the incision area¹⁰.

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Extrusion – Breast-implant extrusion, or breast implant exposure, occurs when the breast skin and tissues that are holding the implant fail, causing the implant to protrude through the skin and become exposed. It happens in less than 2% of patients; it is shortly after breast augmentation or down the road. Breast-implant extrusion can occur for various reasons, such as improper wound healing due to an infection, trauma, too little soft-tissue coverage, oversized implant coupled with too little tissue coverage, or lack of blood supply. Breast-implant extrusion calls for surgery and removal of the implant¹⁴.

Surgical Wound Dehiscence – Surgical wound dehiscence (SWD) is the separation of the margins of a closed surgical incision that has been made in the skin with or without exposure or protrusion of underlying tissue, organs or implants. Separation may occur at single or multiple regions, involve the incision's full length, and affect one or more tissue layers. A dehiscence incision may or may not display

4 Headon H, Kasem A, Mokbel K. Capsular Contracture after Breast Augmentation: An Update for Clinical Practice. *Arch Plast Surg.* 2015 Sep;42(5):532-43. doi: 10.5999/aps.2015.42.5.532. Epub 2015 Sep 15. PMID: 26430623; PMCID: PMC4579193.

5 Song JH, Kim YS, Jung BK, Lee DW, Song SY, Roh TS, Lew DH. Salvage of Infected Breast Implants. *Arch Plast Surg.* 2017 Nov;44(6):516-522. doi: 10.5999/aps.2017.01025. Epub 2017 Oct 27. PMID: 29076316; PMCID: PMC5801791.

6 Franchelli S, Pesce M, Baldelli I, Marchese A, Santi P, De Maria A. Analysis of clinical management of infected breast implants and of factors associated to successful breast pocket salvage in infections occurring after breast reconstruction. *Int J Infect Dis.* 2018 Jun;71:67-72. doi: 10.1016/j.ijid.2018.03.019. Epub 2018 Apr 13. PMID: 29660396.

7 Sood A, Kotamarti VS, Therattil PJ, Lee ES. Sclerotherapy for the Management of Seromas: A Systematic Review. *Eplasty.* 2017 Aug 28;17:e25. PMID: 28890747; PMCID: PMC5575675.

8 Gioacchini M, Bottoni M, Grasseti L, Scalise A, Di Benedetto G. A simple, reliable, and inexpensive method for seroma drainage. *Arch Plast Surg.* 2015 May;42(3):361-2. doi: 10.5999/aps.2015.42.3.361. Epub 2015 May 14. PMID: 26015895; PMCID: PMC4439599.

9 Sforza M, Husein R, Atkinson C, Zaccheddu R. Unraveling Factors Influencing Early Seroma Formation in Breast Augmentation Surgery. *Aesthet Surg J.* 2017 Mar 1;37(3):301-307. doi: 10.1093/asj/sjw196. PMID: 28207027.

10 FDA. Risks and Complications of Breast Implants. Published 2020. <https://www.fda.gov/medical-devices/br-ast-implants/risks-and-complications-breast-implants>

11 Hillard C, Fowler JD, Barta R, Cunningham B. Silicone breast implant rupture: a review. *Gland Surg.* 2017 Apr;6(2):163-168. doi: 10.21037/gls.2016.09.12. PMID: 28497020; PMCID: PMC5409893.

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14 Breast Cancer. Breast Implant Extrusion. Published 2020. [https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after, This is called "extrusion."](https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after, This is called)

clinical signs and symptoms of infection.

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Double Capsule – The double capsule refers to the finding of two distinct capsular layers, separated by an inter-capsular space (ICS), around a breast implant. Although rare, double capsules can occur after breast-implant surgery. This condition's etiopathology is still undefined, but two main hypotheses could explain this complication's development. The first is associated with the periprosthetic capsule's mechanical delamination that creates an intercapsular space (ICS) resulting from fractures after shearing forces are applied between the inner capsule-prosthesis complex and outer capsule. The second hypothesis concerns the development of periprosthetic seromas developing around the capsule. Double capsules may be partial or complete. Clinical signs may vary from asymptomatic to firmness of the implant, discomfort, change in shape or position of the implant, and pain.

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Asymmetry – Preoperative asymmetries include areolas in different positions medially or regarding height, different breast shapes (e.g., one round and the other tuberos), or different breast sizes. These asymmetry types should be differentiated from a postoperative aesthetic difference in the two breasts produced by factors described previously, such as a fall of the fold, a high implant, or a rotation of the implant. Asymmetries caused by the implant's unequal fitting or by creating different sub-mammary grooves. They can be prevented using adequate preoperative planning, correct dissection of the pockets, and comparing the two breasts after fitting the implants. It is possible that after a breast augmentation, small deformities in the wall of the thorax or a morphologic mammary disorder may become much more evident. For this reason, the predicted correction of these anomalies should be discussed with the patient before her operation¹⁷.

Rippling/Wrinkling – Rippling is the cutaneous manifestation, visible or palpable, of the implant ripples and edge that are typically most apparent when the patient bends forward. In situations where the implant's soft tissue coverage is insufficient, these harmful effects become more apparent. Risk factors for rippling are related to the breast-tissue quality and low cohesivity of the implanted gel. Adequate

coverage over the implant ripples is a mandatory element in preventing rippling or implant edge visibility.

Ptois – The "waterfall effect" is a descriptive term to indicate sliding ptosis of parenchymal breast tissue over a fixed or encapsulated implant. It occurs more frequently than surgeons anticipate and especially over the long term after augmentation. Certain breast implants are more prone to contribute to this problem as are implants placed in submuscular pockets that ride high, especially in women with anatomical musculoskeletal variance or asymmetry¹⁸.

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Patients should be advised that their risk of future complications increases with revision surgery compared to primary augmentation or reconstruction surgery. For example, the risk of severe capsular contracture doubles for both augmentation and reconstruction patients with implant replacement compared to first-time implantation. There is a risk of an accidental compromise of implant-shell integrity during reoperation, potentially leading to product failure.

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15 Jong, Justin MD*, Gabriel, Allen MD, FACS†; Trekell, Melissa MD†; Lawser, Amy S, MSN, RNS; Heidel, Eric PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T, MD, FACS** Cohesive Round Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321 doi: 10.1097/GOX.00000000000003321

16 Karan Chopra, MD, Arvind U. Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, W. Grant Stevens, MD. Techniques to Repair Implant Malposition after Breast Augmentation: A Review. Aesthetic Surgery Journal, Volume 36, Issue 6, June 2016, Pages 660–671. <https://doi.org/10.1093/asj/sjv261>

17 Breast Asymmetry (2019). Accessed on March 26, 2021. <https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/asymmetry>

18 Frame J. The waterfall effect in breast augmentation. *Gland Surg.* 2017 Apr;6(2):193-202. doi: 10.21037/igs.2016.10.01. PMID: 28497023; PMCID: PMC5409900

19 Lee Y, Song SE, Yoon ES, Bae JW, Jung SP. Extensive silicone lymphadenopathy after breast implant insertion mimicking malignant lymphadenopathy. *Ann Surg Treat Res.* 2017 Dec;93(6):331-335. doi: 10.4174/ast.2017.93.6.331. Epub 2017 Dec 1. PMID: 29250513; PMCID: PMC5729128.

20 Araco A, Araco E, Sorige R, Gravante G. Sensitivity of the nipple-areola complex and areolar pain following aesthetic breast augmentation in a retrospective series of 1200 patients: Periareolar versus submammary incision. *Plast Reconstr Surg.* 2011;128(4):984-989. doi: 10.1097/PRS.0b013e3182268473

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Interference With Mammography – The importance of mammography should be emphasized. Patients should be instructed to inform their examiners about the presence, type and placement of their implants and request diagnostic mammography rather than screening mammography²¹. Breast implants may complicate the interpretation of mammographic images by obscuring underlying breast tissue and/or compressing overlying tissue. Pre- and post- surgical mammography may be performed to determine a baseline for routine future studies in augmentation patients.

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Interference With Magnetic Resonance Imaging (MRI) – Motiva Implants[®] with a microtransponder are considered MRI conditional. During the MRI study, the microtransponder can create an MRI artifact immediately around the microtransponder (known as an artifact) that can prevent radiologists from being able to see parts of the implant's footprint and parts of the patient's tissue.

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The calculated risk of missing a shell rupture due to the artifact is 1 for every 166,000 units of Motiva Implants[®] with a microtransponder.

The risk of missing breast-cancer detection due to the artifact has been determined to be one high-risk patient with a cancer recurrence for every 596 high-risk patient MRI screening exams performed on patients with Motiva Implants[®] with a microtransponder. When MRI is used in combination with ultrasound (US) to screen the high-risk patient group, it would take 17,892 MRI and US combination screening exams before a patient with cancer recurrence is likely to be missed (false negative).

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23. The NSW Breast Cancer Institute (2006). Breast cysts - An information guide for patients

24. FDA. Risks and Complications of Breast Implants. Published 2020. <https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

25. Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. *Ann Intern Med*. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. Epub 2015 Nov 10. PMID: 26550776.

26. Helmich LR, Lipworth L, McLaughlin JK, Fris S. Breast implant rupture and connective tissue disease: A review of the literature. *Plast Reconstr Surg*. 2007;120(7 SUPPL):1162-69. doi: 10.1097/01.prs.0000286664.50274.f2

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In January 2020, the FDA emitted a medical device report for breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). It included 733 U.S. and global reports of BIA-ALCL of which 496 (68 %) were related to textured breast implants and 28 (4 %) to smooth implants. This data is very similar to the data reported in the previous year, which included a total of 573 BIA-ALCL reports of which 385 (67 %) were related to textured breast implants and 26 (5 %) to smooth implants texture, and one had a history of smooth implant and no known textured implant.²⁹

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28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

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30 Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plastic and reconstructive surgery.* 2017;139(5):1042-1050. doi: 10.1097/PRS.00000000000003282

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31 Cordeiro PG, Ghione P, NIA A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery* : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery.* 2017;140(4):645-654. doi:10.1097/PRS.00000000000003654

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
Breast implants may be positioned in several different tissue planes, described as the pocket location. The selection of this pocket is often a key process in obtaining the desired outcome in individual patients. A pocket that matches the specific conditions of the patient's physical/historical characteristics while maintaining a low-profile risk of complications is crucial. Surgeons should consider the use of published principles for implant selection that minimize the risks of both short and long-term complications.

Incision Site Selection

- The periareolar incision is usually more concealed but may considerably reduce the possibility of future breastfeeding compared to other incision sites. A periareolar incision may be associated with a higher risk of changes in nipple sensation.
- The inframammary incision is generally less concealed than the periareolar, but it is associated with fewer breastfeeding difficulties.
- The axillary incision is the least concealed of all incision sites.

Subpectoral

(under the chest muscle)

PECTORALIS MUSCLE

Subglandular

(under your mammary tissue but over the fascia layer)

FASCIA MAMMARY GLAND

Subfascial

(under both the mammary tissue and the fascia layer)

PECTORALIS MUSCLE FASCIA

Dual Plane

(partly under both the pectoral muscle and mammary gland)

PECTORALIS MUSCLE

Implant Placement Selection

Submuscular/Subpectorol—Benefits of this technique have been reported to include, but are not limited to, greater soft tissue coverage, lower risk of implant visibility, palpability, capsular contracture and infection, less interference with mammography.

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Subglandular—The breast implant is placed behind the breast parenchyma (Subglandular).

Indications for subglandular or subfascial pocket locations may include breast implant primary augmentation, revision augmentation/reconstruction surgeries, patients requiring a pocket exchange from a previous dual plane, submuscular plane, or subfascial plane after revision surgery for a given complication.

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Dual plane—This technique combines varying degrees of soft tissue coverage between the subpectorol placement and the subglandular placement; this enables the surgeon to adjust to the needs of tissue coverage for each patient. The Dual Plane technique is divided into III types. Dual plane I, II, and III.

Indications for dual-plane pocket location may include breast implant primary augmentation, revision augmentation/reconstruction, and primary reconstruction (single/two-stage) surgeries.

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- A backup implant should also be available.
- Be advised that excessive force should not be applied to a small area of the shell during device insertion. Instead, apply force over an area of the implant as large as possible.
- The incision should be of appropriate length to accommodate the volume and profile of the implant with a highly cohesive gel, reducing potential, excessive stress to the implant when inserting it. Forcing implants through a very small opening may damage the implant's gel and cause possible ruptures or gel fractures. In the event of gel fracture during implantation, do not insert the implant into the patient's body and replace it with a new one instead.
- Insufficient pocket dissection increases the risk of rupture and incorrect positioning of the implant. A well-defined, dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.
- All devices should be thoroughly assessed for gel fracture, gel bubbles or any other device failures before being inserted into the patient's body during breast-implant surgery. Do not insert an implant with gel fractures at any point; instead, replace the implant with a new one.
- During explantation, surgeons must assess the breast implant's integrity intraoperatively to identify the presence or absence of rupture, gel fracture and gel migration. In the case of a device failure, the implant should be returned to Establishment Labs for revision.
- **DO NOT** use lubricants during placement since they increase the risk of pocket contamination and may affect the tissue-capsule interface.

- **DO NOT** damage the breast implant with sharp surgical instruments, such as needles and scalpels, blunt instruments (such as clamps and forceps), or over-handling and manipulation during introduction into the surgical pocket.
- **DO NOT** use excessive force during breast-implant placement.
- **DO NOT** manipulate the implant for radial expansion, compression or dissection of the pocket.
- **DO NOT** use more than one implant per breast pocket.

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Establishment Labs strongly recommends that all patients receiving silicone gel-filled breast implants participate in the Establishment Labs device-tracking program, entering their implant information at motiva_health/motivaimagine/. This will help ensure that Establishment Labs has a record of each patient's contact information to be contacted in the event of a field action or other situations with the implants of which they should be made aware.

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³⁵ Establishment Labs. Internal document. RTL-001020 Rationale for Motiva Implants Lifespan.

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The complete terms, conditions and limitations of the Establishment Labs Always Confident Warranty® can be reviewed on the website motiva.health or can be provided by the local Establishment Labs representative. No warranty or program of Establishment Labs covers any costs, fees or expenses related to any medical treatment and/or the surgical replacement of the implants. Establishment Labs shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. In the event that Establishment Labs determines the product was

defective when shipped by Establishment Labs, it shall be the sole responsibility of Establishment Labs to replace the product. Establishment Labs S.A. assumes no further liability. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law, or otherwise, including, but not limited to, any implied warranties of merchantability, suitable for use or performance.

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EUROPEAN REPRESENTATIVE

Emergo Europe:

Westervoortseidijk 60, 6827 AT Arnhem, The Netherlands

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

28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/as/sjz19729 E29 Breast Implants: Evidence Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Último acceso: 28 de febrero de 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breast-implant-associated-anaplastic-large-cell-lymphoma>

21 Sá Dos Reis C, Gremion I, Richli Meystre N. Study of breast implants mammography examinations for identification of suitable image quality criteria. *Insights Imaging*. 3 ene. 2020 ;11(1):3. doi: 10.1186/s13244-019-016-5. PMID: 31900684; PMCID: PMC6942083.

22 Hölmling LR, Mellemkjaer L, Gunnarsdóttir KA, Tange UB, Krag C, Møller S, McLaughlin JK, Olsen JH. Stage of breast cancer at diagnosis among women with cosmetic breast implants. *Br J Cancer*. 24 mar 2003;88(6):832-8. doi: 10.1038/sj.bjc.6600819. PMID: 12644818; PMCID: PMC2377078.

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31 Cordeiro PG, Ghione P, Nj A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. Journal of plastic, reconstructive & aesthetic surgery : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. Annals of surgery. 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

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WARNING: Breast implants have been associated with the development of cancer of the immune system called Breast Implant Associated Anaplastic Large Cell Lymphoma (BIA-ALCL). More detailed information can be found inside on BIA-ALCL section.

Patients with breast implants should have regular clinical follow-up with their surgeon.

Directions For Use

Sterile Silicone Breast Implants

Motiva Implant Matrix®

Establishment Labs

INTRODUCTION

The aim of this Product Insert is to provide an overview of essential information about Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix®, including device description, indications for use, contraindications, warnings, precautions, relevant topics that must be discussed with the patient, adverse events, other reported conditions, returned goods policy, product evaluation, warranty and medical device reporting.


DEVICE DESCRIPTION

Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® are mammary augmentation/ reconstruction devices constructed of successive cross-linked layers of silicone elastomer and a low diffusion barrier shell technology that provide these implants their elasticity and integrity. All implants are composed of the above described shell, a patch, and silicone gel fill. The shell is filled with ProgressiveGel™ PLUS or ProgressiveGel™ ULTIMA®, an Establishment Labs silicone gel proprietary formula. All silicone raw materials are supplied by a FDA approved source in the USA.

The following are the reference ranges of Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix®:

Motiva Implant Matrix® - Round								
Base (cm)	MINI		DEMI		FULL		CORSE	
	P(cm)	V (cc)	P(cm)	V (cc)	P (cm)	V (cc)	P (cm)	V (cc)
8.5	2.2	105	3.1	135	3.5	145	4.0	180
9	2.3	125	3.3	155	3.7	175	4.2	210
9.5	2.4	140	3.4	180	3.9	205	4.5	240
9.75	2.4	150	3.4	190	4.0	220	4.6	260
10	2.5	160	3.5	205	4.1	235	4.8	280
10.25	2.5	170	3.5	215	4.2	255	4.9	300
10.5	2.6	185	3.6	230	4.3	275	5.1	325
10.75	2.6	205	3.7	245	4.4	295	5.2	350
11	2.7	220	3.8	265	4.5	315	5.4	380
11.25	2.7	230	3.8	285	4.6	335	5.5	410
11.5	2.8	245	3.9	300	4.7	355	5.7	440
11.75	2.8	260	3.9	320	4.8	375	5.8	475
12	2.9	275	4.0	340	4.9	400	6.0	510
12.25	2.9	290	4.0	360	5.0	425	6.1	550
12.5	3.0	310	4.1	380	5.1	450	6.3	590
13	3.1	340	4.3	425	5.3	500	6.6	650
13.5	3.2	400	4.4	475	5.5	550	6.9	725
14	3.3	430	4.5	525	5.7	625	7.2	825
14.5	3.4	475	4.6	575	5.9	700	7.5	925
15	3.5	525	4.8	625	6.1	775	7.8	1050

P=Projection
V=Volume

INDICATIONS

Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® are indicated for the following procedures in female patients:

- **Breast augmentation for women of at least 18 years old**, including previous augmentation to increase the breast size and revision surgery to correct or improve the result of a previous breast augmentation surgery.
- **Breast Reconstruction.** Breast reconstruction, including previous reconstruction to replace breast tissue that has been removed due to cancer or trauma, or that has failed to develop properly due to a severe breast anomaly, as well as revision surgery to correct or improve the results of a previous breast reconstruction surgery.

CONTRAINDICATIONS

Breast augmentation with implants is contraindicated in:

- Women with existing carcinoma of the breast, without mastectomy.
- Women with advanced fibrocystic disease considered premalignant, without accompanying subcutaneous mastectomy.
- Women with active infections.
- Women who are currently pregnant or nursing.
- Women with any disease, including uncontrolled diabetes, which is clinically known to impact wound-healing ability.
- Women who show tissue characteristics clinically incompatible with mammoplasty, such as tissue damage resulting from radiation, inadequate tissue, compromised vascularity or ulceration.
- Women with any condition – or treatment – determined by the surgeon to constitute an unjustifiable surgical risk (e.g., unstable cardiovascular disease, coagulopathies, chronic pulmonary problems, etc).

WARNINGS

Care during surgical insertion and subsequent procedures:

- Do not allow sharp instruments, such as scalpels or needles, to come in contact with the device during the implantation or other surgical procedures. Patients should be instructed to inform other treating physicians to also observe this warning.
- Do not immerse the implant in iodine solution. If iodine solution is used in the pocket, make sure that it is rinsed thoroughly with deionized water so that no residual solution remains in the pocket.
- Do not allow the implant to come in contact with cauterization devices.
- Do not alter the implant or attempt to repair or insert a damaged implant.
- Ensure no excessive force is applied to a very small area of the shell during insertion of the device throughout the incision. Instead, apply force over as large an area of the implant as possible during insertion. Excessive forces can lead to implant failure by either gel fracture or implant rupture.

- The incision should be of appropriate length to accommodate the volume and profile of the implant. This will reduce the potential for creating excessive stress to the implant when inserting it. Forcing implants through a very small opening may result in local weakening of the breast implant shell potentially leading to shell damage, gel fracture and possible implant rupture.
- Periareolar and axillary incision sites may make the insertion more difficult, increasing the risk of damage to the implant. Periareolar incision may considerably reduce the possibility of future breastfeeding.
- Do not use the periumbilical approach to place the implant.
- Avoid creating wrinkles or folds in the device during the insertion. It is recommended to run a finger around the implant before closing to make sure the implant is flat.
- Do not place more than one implant per breast pocket.
- Do not treat capsular contracture by closed capsulotomy or forceful external compression, which will likely produce implant damage, rupture, folds, and/or hematoma.
- Procedures such as open capsulotomy, breast pocket revision, hematoma/seroma aspiration, biopsy, and lumpectomy might result in damage to the implant shell, so they must be carefully performed. Care should be taken when re-positioning the implant during subsequent procedures to avoid contamination of the implant. Use of excessive force during any subsequent procedure can contribute to local weakening of the breast implant shell, potentially leading to shell damage and possible implant rupture.
- Do not re-use or re-sterilize any product that has been previously implanted. Breast implants are intended for single use only.
- Do not use microwave diathermy in patients with breast implants, as it has been associated with tissue necrosis, skin erosion, and implant extrusion.

PRECAUTIONS

1. Specific Populations

Safety and effectiveness of breast augmentation surgery have not been established for the following populations and/or conditions:

- Patients with autoimmune diseases (e.g., lupus, scleroderma).
- Patients whose immune system is compromised (e.g., currently receiving immunosuppressive therapy such as steroids).
- Patients with conditions or medications that may interfere with wound healing ability (e.g., poorly controlled diabetes, or corticosteroid therapy) or blood clotting (e.g., concomitant Warfarin therapy).
- Patients with reduced blood supply to breast or overlying tissue.
- Patients undergoing radiation therapy.
- Women with ptotic breasts where nipples fall below the inframammary fold, without concurrent mastopexia.
- Previous repeated contour correction failures.
- Patients with clinical diagnosis of depression or other mental health disorders, including BDD (body dysmorphic disorder) and eating disorders. The patient should be advised to discuss any history of mental health disorders with her surgeon prior to surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until stabilization of these conditions prior to undergoing breast implant surgery.

- There may be other patients with complicated medical histories, who are judged to present risk factors that might interfere with the safety and effectiveness of the breast implant surgery. As with any surgical procedure, the patient's medical history should be carefully reviewed to ensure that she is an appropriate candidate for breast implant surgery.

2. Surgical Precautions

Preliminary product examination– Immediately before insertion, examine the device by gently manipulating it while carefully checking for rupture, gel fracture, leakage sites or particulate contamination.

Surgical technique and implant selection– There are several surgical techniques that can be used to perform the implantation of a silicone gel-filled breast implant. Therefore, the surgeon is advised to use his/her clinical judgment in choosing the procedure that is best for the patient, consistent with this product insert. The incision should be of appropriate length to accommodate the volume and profile of the implant with highly cohesive gel. This will reduce the potential for creating excessive stress to the implant when inserting it. Forcing implants through a very small opening may result in damage to the implants gel and possible ruptures or gel fractures. After setting realistic aesthetic goals that assure mutual understanding between doctor and patient, the surgeon must choose from current and accepted surgical techniques to minimize the incidence of adverse reactions and achieve the best results.

The implant size should be consistent with the patient's chest wall dimensions, including base width measurements, characteristics of the tissue and projection of the implant.

Textured implants, larger implants, subglandular placement, and an insufficient amount of tissue available to cover the implant may cause them to be more palpable.

Implants of larger sizes may increase the risk of complications such as extrusion, hematoma, infection, palpable implant folds, and visible skin wrinkling.

CAUTION: Only surgeons with qualified training and certified by the corresponding national medical board of their countries should use this product. The use of this product by unqualified practitioners may result in extremely poor aesthetic outcomes and serious adverse effects.

RELEVANT TOPICS THAT MUST BE DISCUSSED WITH THE PATIENT

Patient Counseling Information

This document and the information for the patient must be carefully reviewed before counseling a patient about Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® and breast augmentation surgery. Doctors must thoroughly read and understand the contents of these document and make sure that any questions or concerns have been resolved prior to proceeding with use of the device. Breast implant surgery is an elective procedure, and the patient must understand its potential risks and benefits in order to make an informed decision. For this reason, the patient should be instructed to read the document called "**DOC-001004 Information for the patient breast augmentation and reconstruction with Motiva Implants®**". The doctor must discuss with the patients sections on warnings, contraindications, precautions, important factors to consider, complications and all other aspects of the document. The physician should inform the patient about the potential complications and that medical management of serious complications may include additional surgery and explanation.

Informed Consent

Each patient should be given the Establishment Labs "**DOC-001004 Information for the patient breast augmentation and reconstruction with Motiva Implants®**"

during her surgical consultation. It is the surgeons' responsibility to ensure this happens and it is a requirement for the use of the device. The patient must be given enough time to read and completely understand the information regarding the risks, benefits and recommendations associated with silicone gel-filled breast implant surgery.

In order to document a successful informed decision process, the patient, a witness and the surgeon should sign the **"Informed Consent Document"**, which will be part of the patient's medical file.

Some of the relevant topics patients need to be aware of when considering the use of silicone gel-filled breast implants are:

Rupture– Breast implants can rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation, but it is more likely to happen when the implant has been in place for a long time. Rupture of a silicone gel-filled breast implant is most often silent (the patient does not experience any symptoms and there are no physical signs of changes with the implant) rather than symptomatic. Therefore, patients should be advised to have regular MRIs over their lifetime to screen for silent rupture even if they are not having any apparent problems. The first MRI should be performed at 3 years postoperatively, then regularly at 2-year intervals, and the images submitted to the treating surgeon. Patients should be provided with a list of radiology centers with experience with breast implant MRI to scan for signs of rupture. The importance of these MRI evaluations should be emphasized. If rupture is noted on an MRI, the patient should be strongly encouraged to have her implant removed.

Gel Fracture– Gel fracture can occur with cohesive silicone and occur most frequently as a consequence of subjecting the implant to excessive compression forces during the implantation. Alternatively, gel fracture can occur due to the development of capsular contracture and may result in device distortion.

Mammography– Patients should be advised to have routine mammography performed according to their surgeon's recommendations. The importance of these exams should be emphasized. Patients should be instructed to inform their examiners about the presence, type, and placement of their implants, and to request a diagnostic mammography, rather than a screening mammography. Breast implants may complicate the interpretation of mammographic images by obscuring underlying breast tissue and/ or by compressing overlying tissue. Accredited mammography centers, technicians with experience in examining patients with breast implants and the use of displacement techniques are needed to adequately visualize breast tissue in the implanted breast. Pre and post-surgical mammographies may be performed to determine a baseline for routine future studies in augmentation patients.

Explantation– Implants are not lifetime devices, and there is a possibility that patients will undergo implant removal(s), with or without replacement, over the course of their life. When implants are explanted without replacement, changes to the patient's breasts may be irreversible. Complication rates are higher following revision surgery (removal with replacement).

Reoperation– Rupture, unacceptable cosmetic outcomes (dimpling, wrinkling, and other potentially permanent cosmetic changes of the breast) and other complications may require additional surgeries to the patient's breasts. Patients should be advised that their risk of future complications increases with revision

surgery as compared to primary augmentation or reconstruction surgery. For example, the risk of severe capsular contracture doubles for both augmentation and reconstruction patients with implant replacement compared to first time implantation. There is a risk of an accidental compromise of implant shell integrity during reoperation, potentially leading to product failure.

Infection– Signs of acute infection reported in association with breast implants include edema, erythema, tenderness, pain, and fever. As with other invasive surgeries, Toxic Shock Syndrome (TSS), a life-threatening condition, has been reported in rare instances following breast implant surgery. Symptoms of TSS occur suddenly and can include high fever (102°F, 38.8°C or higher), vomiting, diarrhea, sunburn-like rash, red eyes, dizziness, lightheadedness, muscle aches and drops in blood pressure, which may cause fainting. Patients should immediately contact their physician for diagnosis and treatment if any of these symptoms occur.

Breast Examination Techniques– Patients should perform breast self-examinations monthly and be shown how to distinguish the implant from their breast tissue. The patient should not manipulate or squeeze the implant excessively. The patient should be told that the presence of lumps, persistent pain, swelling, hardening, or change in the implant shape could suggest symptomatic rupture of the implant. If the patient has any of these signs, she should be advised to report them, and possibly have an MRI evaluation to screen for rupture.

Lactation– Breast implant surgery may interfere with the ability to successfully breast feed, either by reducing or eliminating milk production. Particularly, periareolar incision may considerably reduce the possibility of breast-feeding.

Avoiding Damage During Treatment– Patients should inform other treating physicians of the presence of implants to minimize the risk of damage to the implants.

Topical Medications– The patient should consult a physician or a pharmacist before the use of topical medicines (e.g. steroids) in the breast area.

Trauma– The patient should consult the surgeon or physician if any complications are suspected, in particular in the case of trauma or compression caused, for example, by extreme massaging of the breast region, by some sport activities or by using seat belts.

Smoking– Smoking may interfere with the healing process.

Radiation to the Breast– Establishment Labs has not tested the in vivo effects of radiation therapy in patients who have breast implants. The literature suggests that radiation therapy may increase the likelihood of capsular contracture, necrosis, and implant extrusion.

Insurance Coverage– Patients should check with their insurance company regarding coverage issues before undergoing surgery.

Mental Health and Elective Surgery– It is important that all patients seeking an elective procedure such as breast augmentation have realistic expectations that focus on improvement rather than perfection. Ask the patient to openly discuss, prior to surgery, any history of depression or other mental health disorders.

Postoperative Care:

The patient should be advised that she will likely feel tired and sore for several days following the operation, and that her breasts may remain swollen and sensitive to physical contact for a month or longer. She also could experience a feeling of tightness in the breast area as the skin adjusts to the new breast size. The patient should avoid any strenuous activities for at least a couple of weeks but should be able to return to work within a few days. Breast massage may also be recommended as appropriate.

Life expectancy of the breast implant:

Breast implants are not lifetime devices. However, the life expectancy of a silicone breast implant cannot be precisely estimated, as there are many factors beyond the manufacturer's control that can affect the longevity of a device. The time period varies from woman to woman. Some women could need replacement surgery few years after the augmentation procedure while others can have their implants intact for 10-20 years. Therefore, the life expectancy of the implant cannot be guaranteed.

ADVERSE EVENTS

If any of the following or other adverse events occurs fill out a complaint notification form, providing all available information on patients, product information, a reason for the complaint and a summary of the event and send it to <https://motiva.health/surgeons/>.

Because breast implant surgery is more often performed using general anesthesia, it is associated with the same risks as other invasive surgical procedures. After breast implant surgery, patients might experience swelling, hardness, discomfort, itching, allergies, bruising, twinges and pain over the first few weeks. Potential adverse events that may occur with silicone gel-filled breast implant surgery include:

Capsular Contracture

Normally, capsules of collagen fibers form as an immune response around a foreign body, such as a breast implant, tending to isolate it. Capsular contracture occurs when the capsule tightens and squeezes the implant. This can cause the implant to turn rigid (from slightly firm to quite hard) and the firmest ones can cause varying degrees of discomfort, pain and palpability. In addition to the firmness, capsular contracture can result in a deformed breast, visible surface wrinkling and/or displacement of the implant. Detection of breast cancer by mammography may also be more difficult. Capsular contracture may be more common following infection, hematoma, and seroma, and the chance of it happening may increase over time. Capsular contracture occurs more commonly in patients undergoing revision surgery than in patients undergoing primary implantation surgery. Capsular contracture is a risk factor for implant rupture, and it is the most common reason for reoperation in augmentation and reconstruction patients.

Capsular contracture is graded into 4 levels depending on its severity. Baker Grade I: the breast is normally soft and looks natural; Baker Grade II: the breast is a little firm but looks normal; Baker Grade III: the breast is firm and looks abnormal; Baker Grade IV: the breast is hard, painful, and looks abnormal. Patients should also be advised that additional surgery might be needed in cases where pain and/or firmness are severe (Baker Grades III or IV) and that capsular contracture may happen again after additional surgeries.

Correction of capsular contracture may require surgical removal or release of the capsule, or removal and possible replacement of the implant itself. Closed capsulotomy (external manipulation of the capsule in order to "pop" the tissue

capsule and open it up) used to be a common procedure for treating capsular contracture, but most manufacturers, including Establishment Labs, contraindicate it because it can cause implant rupture.

Rupture

Breast implants can potentially remain intact for decades in the body, but all such devices will fail at some point.

Breast implants rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation, but it is more likely to occur the longer the implant is in place. The following may cause implants to rupture: damage by surgical instruments, implant stress and weakening during implantation, age and design of the implant, submuscular rather than subglandular location, occurrence of post-operative hematomas or seromas, folding or wrinkling of the implant shell, excessive force to the chest (e.g., during closed capsulotomy, which is contraindicated), trauma, compression during mammographic imaging, and severe capsular contracture.

Silicone gel-filled implant ruptures are most often silent. (MRI examination is currently the best screen method for silent rupture.) This means that most of the time neither doctor nor patient will know if the implant has a tear or hole in the shell. This is why a first MRI is recommended after 3 years, and then at regular intervals every 2 years thereafter, to screen for ruptures. Sometimes there are symptoms associated with gel implant rupture, such as lumps surrounding the implant or in the armpit, change or loss of size or shape of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast.

When MRI findings of rupture are found, or if there are signs or symptoms of rupture, the implant should be removed, with or without replacing it. If the patient develops symptoms that suggest implant rupture, she should be told to have an MRI evaluation to screen for rupture.

If rupture occurs, silicone gel may either remain within the scar tissue capsule surrounding the implant (intracapsular rupture), move outside the capsule (extracapsular rupture), or move beyond the breast (migrated gel). There is also a possibility that rupture may progress from intracapsular to extracapsular and beyond.

Below is a summary of information related to the health consequences of implant rupture, which **have not been fully established**, in women who had a variety of implant models from different manufacturers.

- Local breast complications that were associated with rupture in the literature include breast firmness, a change in breast shape or size, and breast pain. These symptoms are not specific to rupture and may be also experienced by women who have capsular contracture.
- There have been rare reports of gel migration to nearby tissues such as the chest wall, armpit, or abdominal wall, and to more distant locations down the arm or into the groin. This has led to nerve damage, granuloma formation and/or breakdown of tissues in direct contact with the gel in a few cases. There have been reports of silicone presence in the liver of patients with silicone breast implants. Movement of silicone gel material to lymph nodes in the armpit also has been reported, even in women without evidence of rupture, leading to lymphadenopathy.
- Concerns have been raised over whether ruptured implants are associated with the development of connective tissue or rheumatic diseases and/ or symptoms such as fatigue and fibromyalgia. A number of epidemiology studies have evaluated large populations of women with breast implants from a variety of manufacturers and implant models. These studies do not support an association of breast implants and rheumatic disease.

Gel Fracture

Gel fracture is defined as a fissure or crack in the gel that occurs when the filler is forcibly separated by exceeded intrinsic forces. As a result, the shape is irrevocably lost requiring replacement. (Jill Baker et al,2016).

Gel fractures can occur with cohesive silicone and occur most frequently as a consequence of subjecting the implant to excessive compression forces during the implantation. Alternatively, gel fracture can occur due to the development of capsular contracture and may result in device distortion.

Be advised that during device insertion excessive force should not be applied to a small area of the shell. Instead, apply force over as large an area of the implant as possible.

The incision should be of appropriate length to accommodate the volume and profile of the implant with highly cohesive gel. This will reduce the potential for creating excessive stress to the implant when inserting it. Forcing implants through a very small opening may result in damage to the implants gel and possible ruptures or gel fractures.

Gel fracture can be detected by ultrasound or magnetic resonance imaging (MRI). The majority of gel fractures are clinically undetectable.

Pain

Most women undergoing augmentation or reconstruction with a mammary implant will experience some post-operative breast and/or chest pain. While this pain normally recedes in most women as they heal after surgery, it can become a chronic problem in other women.

Hematoma, migration, infection, implants that are too large, or capsular contracture can cause chronic pain. Sudden, severe pain may be associated with implant rupture. The surgeon should instruct the patient to immediately report if there is significant pain or if pain persists.

Changes in Nipple and Breast Sensation

Breast surgery can result in an increased/decreased breast and/or nipple sensitivity. Typically the sensation is lost after complete mastectomy where the nipple itself is removed, and can be severely lessened by partial mastectomy. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. While some of these changes can be temporary, they can also be permanent, and may affect the patient's sexual response or ability to nurse.

Infection

Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are harder to treat than infections in tissue without an implant. If an infection does not respond to antibiotics, the implant may have to be removed, and another implant may be placed after the infection is resolved. As with other surgical procedures, toxic shock syndrome in rare instances has been noted in women after breast implant surgery. This is a life-threatening condition and its symptoms include sudden fever, vomiting, diarrhea, fainting, dizziness, and/or sunburn-like rash. Patients should be instructed to contact their doctor immediately for diagnosis and treatment if they have these symptoms.

Hematoma/Seroma

Hematoma is a collection of blood within the space around the implant, and a seroma is a build-up of fluid around the implant. Having a hematoma and/ or seroma following surgery may result in infection and/ or capsular contracture later on. Symptoms from a hematoma or seroma may include swelling, pain, and bruising. If a hematoma or seroma occurs, it will usually be soon after surgery. However, they can also occur at any time after injury to the breast. While the body absorbs small hematomas and seromas, some will require surgery, typically involving draining, and potentially placing a surgical drain in the wound temporarily for proper healing. A small scar can result from surgical draining. Implant rupture also can occur from surgical draining if there is damage to the implant during the procedure.

Breast-feeding

Although most women with breast implants who attempt nursing have successfully breastfed their babies, it is not known if there are increased risks for a woman with breast implants or if the children of women with breast implants are more likely to have health problems. At this time, it is not known if it is possible for a small amount of silicone to pass from the breast implant silicone shell into breast milk during breastfeeding, or what the potential consequences might be.

A periareolar surgical approach may further increase the chance of breastfeeding difficulties. However, the American Academy of Pediatrics has stated that there is no reason why a woman with implants should refrain from nursing.

Calcification

Calcium deposits can form in scar tissue surrounding the implant and may cause pain and firmness, and be visible on a mammography. These deposits must be identified as different from calcium deposits that are a sign of breast cancer. Additional surgery may be necessary to remove and examine calcifications. Calcium deposits also occur in women who undergo breast reduction procedures, in patients who have had hematoma formation, and even in the breasts of women who have not undergone any breast surgery. The occurrence of calcium deposits significantly increases with age.

Delayed Wound Healing

Some patients may experience a prolonged wound healing time. Smoking may interfere with the healing process. Delayed wound healing may increase the risk of infection, extrusion, and necrosis. Wound healing times may vary depending on the type of surgery or incision.

Implant Extrusion

Lack of adequate tissue coverage, local trauma or infection may result in exposure and extrusion of the implant. This has been reported with the use of steroid drugs or after radiation therapy of breast tissue. If tissue breakdown occurs and the implant becomes exposed, implant removal may be necessary, which may result in additional scarring and/ or loss of breast tissue.

Necrosis

Necrosis is the formation of dead tissue around the implant. This may prevent wound healing and require surgical correction and/ or implant removal. Permanent scar deformity may occur following necrosis. Factors associated with necrosis include infection, use of steroids in the surgical pocket, smoking, chemotherapy/ radiation, and excessive heat or cold therapy.

Granulomas

These are benign lumps that can form when body cells surround foreign material such as silicone. Like any lump, it should be further evaluated to rule out a malignancy.

Breast Tissue Atrophy/Chest Wall Deformity

The pressure of the breast implant may cause the breast tissue to thin and shrink (with increased implant visibility and palpability), potentially leading to chest wall deformity. This can occur while implants are still in place or following implant removal without replacement. Either of these conditions may result in additional surgeries and/or unacceptable dimpling/creasing of the breast.

Lymphadenopathy

Literature reports associate lymphadenopathy with both intact and ruptured silicone breast implants. One study reported that armpit lymph nodes from women with both intact and ruptured silicone gel implants had abnormal tissue reactions, granulomas, and the presence of silicone.

These reports occurred in cases of women who had implants from a variety of manufacturers and implant models.

Unsatisfactory Results

Unsatisfactory results such as wrinkling, asymmetry, implant displacement/migration, incorrect size, implant palpability/visibility, scar deformity, and/or hypertrophic scarring, may occur. Some of these results may cause discomfort. Pre-existing asymmetry may not be entirely correctable by implant surgery. Revision surgery could be indicated to increase patient satisfaction, but this involves additional considerations and risks. Careful preoperative planning and surgical technique can minimize but not always prevent unsatisfactory results.

Other Reported Conditions

There have been reports in the literature of other conditions in women with silicone breast implants. Many of these conditions have been studied to evaluate their potential association with breast implants. However, no causal relationship has been established between breast implants and the conditions listed below.

Connective Tissue Disease (CTD)

Since the early 1990s, nearly a dozen comprehensive systemic reviews have been commissioned by government health ministries in several countries to examine the alleged links between silicone gel breast implants and systemic diseases. A clear consensus has emerged from these independent scientific reviews that there is no clear evidence of a causal link between the implantation of silicone breast implants and connective tissue disease.

Cancer

Breast cancer reports in the medical literature reveal that patients with breast implants are not at a greater risk than those without breast implants for developing breast cancer. Some reports have suggested that breast implants may interfere with or delay breast cancer detection by mammography and/ or biopsy; however, other reports in the medical literature indicate that breast implants do not significantly delay breast cancer detection or adversely affect cancer survival prognosis in implanted women. Some studies even suggest lower rates of breast cancer in women with breast implants.

Neurological Disease, Signs, and Symptoms

Some women with breast implants have experienced neurological disturbances (e.g., visual symptoms or alterations in sensation, muscle strength, walking, balance, thinking or memory) or diseases (e.g., multiple sclerosis) and they believe those symptoms are related to their implants. However, there is no

evidence in published literature of a causal relationship between breast implants and neurological disease.

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Small quantities of silicone may diffuse through the elastomer envelope of silicone gel-filled implants. The detection of small quantities of silicone in the periprosthetic capsule, axillary lymph nodes and other distal regions in patients with apparently intact gel-filled implants has been reported in the literature. Some studies on long-term implants have suggested that gel-bleed may contribute to the development of capsular contracture and lymphadenopathy. On the other hand, evidence against gel-bleed being a significant contributing factor to capsular contracture and other local complications is provided by the fact that there are similar or lower complication rates for silicone gel-filled breast implants than for saline-filled breast implants.

Malposition

Malposition of a breast implant is defined as an incorrect placement during surgery or the shifting of the implant from its original position. Malposition has been a frequent event reported due to its multifactorial causes and it can be expected during the lifetime of the device.

The shifting of the implants can be produced by trauma, capsular contracture, gravity, or initial wrong placement. The surgeon must plan the operation carefully and conduct the surgery with a technique that can minimize, but not completely evade the risk of malposition. The risk associated with this event is dissatisfaction with aesthetic outcomes.

The clinical symptoms manifested by the patients are change in breast shape, displacement or sensation of firmness. Revision surgery may be indicated to achieve patient satisfaction. New considerations and risks need to be taken into account before performing a revision surgery.

Bottoming-out

The inferior displacement of a breast implant, increasing the distance between the nipple-areolar complex and the inframammary fold, after breast implant surgery. Risk factors reported in the literature are, but not limited to the quality of the preexisting breast tissue (thin subcutaneous tissue, defective dermal elements, breast tuberosity), the breast implant selection (larger implants), the IMF dissection and the placement of the implant during surgery (submuscular and subglandular planes). The clinical symptoms resulting from a bottoming out implant are asymmetry, upward-pointing nipples, sagging breast, palpable implant, among others. Prevention of this complication is based on the anticipation of the possible causes, for example: a careful and individual assessment of the mammary soft-tissues, a careful implant selection, preparation with a technique that can minimize the risk, and to provide adequate breast support after the surgery. The treatments may vary depending on the severity of the complication and go from a simple sub mammary fixation to the use of additional supporting materials.

Flipping

Anterior/posterior malposition, also called flipping, has been described to occur more frequently with cohesive gel implants. The shape of the breast is lost because the flat base of the implant is positioned anteriorly, deforming the breast of the patient. It has been reported in the literature that the interaction between breast envelopes, physical characteristics of the implant, and the pocket dissection is the cause of malposition. Other theories include the involution of the breast tissue. Regarding the implant characteristics, it has been associated with the presence or absence of texturing, the shape/profile of the implant, and the gel-filling ratio. Other factors such as infection, hematoma, capsular contracture, dissection,

surgeon's experience, physical activity, and external manipulation of the implant could potentially contribute to the development of this complication. The diagnosis is based on clinical evidence: the patient complaints of loss of breast shape, which is confirmed by the inspection of the affected area. MRI or CT imaging to validate the diagnosis can be useful but are not necessary. Flipping can be treated with bimanual manipulation in the office and can be repeated in recurrent cases. However, in some cases, it may be necessary a revision surgery to reduce pocket dimensions.

Interference with Mammography

Breast implants (especially in subglandular placement) may complicate the interpretation of mammographic images by obscuring the underlying breast tissue and/or by compressing overlying tissue. Despite of the fact that the presence of breast implants lessens tissue compression range during mammography, a number of studies looking at breast cancers in women with implants have found no significant difference in stage of disease at time of diagnosis, and prognosis appears to be similar in implanted and nonimplanted patients. Accredited mammography centers, technicians with experience in imaging patients with breast implants, and use of displacement techniques are needed to adequately visualize breast tissue in the implanted breast. Anterior breast tissue is best visualized with displacement views and posterior breast tissue with compression views. The decrease in visible area of 35% with compression views is improved to 25% with displacement views.

Interference with Magnetic Resonance Imaging (MRI)

Sterile Silicone Breast Implants Motiva Implant Matrix® with a *microtransponder* are considered MRI conditional. The *microtransponder*, during an MRI study, can create an MRI image immediately around the *microtransponder* (known as an artifact) that can prevent radiologists from being able to see parts of the implant's footprint and parts of the patient's tissue.

Therefore, there are potential added MRI risks associated with this artifact including, but not limited to, an inadequate evaluation of the implant shell for the detection of rupture or missing a diagnosis of cancer should it obscure a cancer in the artifact area.

Calculated risk of missing a shell rupture due to the artifact is 1 for every 166,000 units of Motiva Implants® with *microtransponder*.

Risk of missing breast cancer detection due to the artifact has been determined to be 1 high-risk patient with a cancer recurrence for every 596 high-risk patient MRI screening exams performed on patients with Motiva Implants® with *microtransponder*. When MRI is used in combination with US, to screen the high-risk patient group, it would take 17,892 MRI and US combination screening exams before a patient with cancer recurrence is likely to be missed (false negative).

Reduction of these risks can be achieved by performing an Ultrasound (US) in addition to the MRI thereby allowing the radiologist to see the area within the artifact produced in the MRI. As such, the radiologist should be informed of the presence of the *microtransponder* and that it is embedded near the patch area inside the breast implant. The presence of an MRI imaging void artifact should be anticipated along with its expected size.

Breast Implant Associated-Anaplastic Large Cell Lymphoma (BIA-ALCL)

BIA-ALCL is a rare type of T-cell lymphoma involving cells of the immune system. The World Health Organization in the year 2016 recognized it as a breast implant-associated disease. The exact number of cases remains difficult to determine due to significant limitations in worldwide reporting and lack of global implant sales data. It has been reported that most data suggest that BIA-ALCL occurs more frequently following implantation of breast implants with textured surfaces rather than those with smooth surfaces.

The French National Agency for Medicines and Health Products Safety (ANSM) has requested manufacturers of textured breast implants to perform biocompatibility testing. Establishment Labs has complied with this request.

There is a significant body of medical literature relating to breast implants and the risk of developing ALCL. According to the FDA, all of the information reviewed as of the date of FDA's March 2017 notice, suggests that "women with breast implants have a very low but increased risk of developing ALCL compared to women who do not have breast implants". Most cases of breast implant associated ALCL are treated by removal of the implant and the capsule surrounding the implant and some cases have been treated by chemotherapy and radiation.

The following are considerations from the FDA to the investigators regarding BIA-ALCL:

If you have patients with breast implants, you should continue to provide them routine care and support including regular clinical follow-up. BIA-ALCL is a very rare condition; when it occurs, it has been identified most frequently in patients undergoing implant revision operations for late onset persistent seroma. Because it has generally only been identified in patients with late onset of symptoms such as pain, lumps, swelling, or asymmetry, prophylactic breast implant removal in patients without symptoms or other abnormality is not recommended.

Current recommendations include the steps below:

Be aware that most confirmed cases of BIA-ALCL have occurred in women with textured breast implants. Provide the manufacturers labeling as well as any other educational materials to your patients before surgery and discuss with them the benefits and risks of the different types of implants.

Consider the possibility of BIA-ALCL when you have a patient with late onset, persistent peri-implant seroma. In some cases, patients presented with capsular contracture or masses adjacent to the breast implant. If you have a patient with suspected BIA-ALCL, refer the patient to an appropriate specialist for evaluation.

When testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out BIA-ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell block immunohistochemistry testing for cluster of differentiation (CD) and Anaplastic Lymphoma Kinase (ALK) markers.

Develop an individualized treatment plan in coordination with the patient's multi-disciplinary care team. Consider current clinical practice guidelines, such as those from the Plastic Surgery Foundation or the National Comprehensive Cancer Network (NCCN) when choosing your treatment approach.

For the latest statistical data on reported cases refer to:

<https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/BreastImplants/ucm239995.htm>

INSTRUCTIONS FOR USE

Single Use

This product is intended to be used only in one patient for a single procedure. DO NOT reuse explanted implants. To reuse a single-use device could expose patients and staff to risks which outweigh the perceived benefits of using such devices. This product is not intended to be reprocessed in any way and/or used again, not even on the same patient. The reuse of single-use devices can affect their safety, performance and effectiveness, exposing patients and staff to unnecessary risk, such as infection, the inability to assure the proper cleaning and decontamination, the presence of residues of cleaning agents, reaction to endotoxins, exposure to other bio-hazards and/ or device failure. This practice may also have legal implications that vary according to each jurisdiction.

Product Traceability

The product traceability stickers, provided with each device and located within the internal product packaging, provide product specific information and should be attached to the patient's chart for identification purposes. Stickers are also available for the Patient ID Card and the hospital files, if applicable. The surgeon should encourage the patient to participate in the Establishment Labs device-tracking program, entering their implant(s) information at www.motivaimplants.com.

This will help ensure that Establishment Labs has a record of each patient's contact information so that they can be contacted in the event of a field action or other problems with the implants of which they should be made aware.

Sterile Product

Product is sterilized by the manufacturer using Dry Heat Sterilization method, each sterile silicone breast implant is supplied in a sealed, double sterile barrier primary package. Use standard procedures to maintain sterility during transfer of the breast implant to the sterile field. Remove the breast implant from their packages in an aseptic environment, using talc-free gloved hands.

Sterility of the implant is maintained only if the thermoform packages, including the package seals, are intact.

DO NOT use the product if the thermoform packages or seals have been damaged.

DO NOT re-sterilize the product.

Avoid prolonged exposure to extreme storage conditions. We recommend keeping these devices at room temperature, at atmospheric pressure, in dry conditions and away from direct sunlight.

DO NOT implant any device that may appear to have particulate contamination, damage, or loss of shell integrity. A sterile back-up implant must be readily available at the time of surgery.

DO NOT implant any device that may appear to have leaks or scratches.

How to Open Sterile Product Package


DO NOT expose the breast implant to talc, sponges, towels, or other contaminants.

1. A non-sterile team member should open the outer package.
2. Remove the inner package and invert it over the sterile field, allowing the sealed inner thermoform package to gently slide into the field.
3. Use the pull-tab to open the lid of the inner thermoform package.
4. Retrieve the breast implant and examine it for any particulate contamination, damage, or loss of shell integrity. If satisfactory, return the breast implant to the inner thermoform tray. At this point, you might slightly rinse the implant with a small amount of saline to remove the static and cover the tray with the lid until implantation to prevent contact with airborne and surgical field particulate contaminants.

Surgical Technique and Implant Selection


There are several surgical techniques that can be used to perform the implantation of silicone gel-filled breast implants. Therefore, the surgeon is advised to use his/her clinical judgment in choosing the procedure that is best for the patient, consistent with this product insert. After setting realistic aesthetic goals that assure mutual understanding between doctor and patient, the surgeon must choose from current and accepted surgical techniques to minimize the incidence of adverse reactions and achieve the best results.

The surgeon must carefully assess the implant size, projection and surface, as well as the incision placement, pocket dissection and implant placement criteria according to the patients' anatomy and desired aesthetic outcomes.

Incision Site Selection

- The periareolar incision is usually more concealed but may considerably reduce the possibility of future breastfeeding as compared to other incision sites. A periareolar incision may be associated with a higher risk of changes in nipple sensation.
- The inframammary incision is generally less concealed than the periareolar, but it is associated with less breast-feeding difficulties.
- The axillary incision is the least concealed of all incision sites.
- The periumbilical approach should not be used with Sterile Silicone Breast Implants Motiva Implant Matrix[®] for a number of reasons, including potential damage to the implant shell.

Implant Placement Selection

- The possible benefits of submuscular placement are that it may result in less palpable implants, less likelihood of capsular contracture, and easier mammographies. This placement may be preferable if the patient has thin or weakened breast tissue. However, submuscular placement is associated with a longer surgical procedure, a more prolonged recovery period and more pain. Also, it may make it more difficult to perform some reoperation procedures.
- Subglandular placement may make surgery and recovery shorter, be less painful, and provide easier access for reoperation than the submuscular placement. However, this placement may result in more palpable implants, greater risk of capsular contracture, ptosis and increased difficulty in imaging the breast with mammography.

- Dual plane placement has been associated by some authors with the benefits of submuscular placement with the advantages of a faster recovery and less pain and postoperative discomfort.

During the Surgical Procedure:

- It is advisable to have more than one size of breast implant in the operating room at the time of surgery to allow for flexibility in determining the appropriate size to be used.
- A backup implant should also be available.
- Be advised that during device insertion excessive force should not be applied to a small area of the shell. Instead, apply force over an area of the implant as large as possible.
- The incision should be of appropriate length to accommodate the volume and profile of the implant with highly cohesive gel. This will reduce the potential for creating excessive stress to the implant when inserting it. Forcing implants through a very small opening may result in damage to the implants gel and possible ruptures or gel fractures. In the event of gel fracture during implantation, do not insert the implant in to the patient's body and replace it with a new one
- Insufficient pocket dissection increases the risk of rupture and incorrect positioning of the implant. A well-defined, dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.
- During breast implant surgery, all devices should be thoroughly assessed for gel fracture, gel bubbles or any other device failures before being inserted into the patient's body. Do not insert an implant with gel fractures at any point, instead replace the implant with a new one.
- During explantation, surgeons must assess intraoperatively the integrity of the breast implant to identify presence or absence of rupture, gel fracture and gel migration. In the case of a device failure the implant should be returned to Establishment Labs for revision.
- **DO NOT** use lubricants during placement since they increase the risk of pocket contamination and may also affect the tissue-capsule interface.
- **DO NOT** damage the breast implant with sharp surgical instruments such as needles and scalpels, blunt instruments such as clamps and forceps, or by over-handling and manipulation during introduction into the surgical pocket.
- **DO NOT** use excessive force during breast implant placement.
- **DO NOT** manipulate the implant for radial expansion, compression or dissection of the pocket.
- **DO NOT** use more than one implant per breast pocket.

Maintaining Hemostasis/Avoiding Fluid Accumulation

The risk of postoperative hematoma and seroma may be reduced by carefully The risk of postoperative hematoma and seroma may be reduced by carefully handling the hemostasis during surgery, and possibly also by postoperative use of a closed drainage system. Persistent or excessive bleeding must be controlled before implantation.

Any postoperative evacuation of hematoma or seroma must be conducted with care to avoid contamination or damage to the breast implant.

Instructions and precautions for removal

Among the most common reasons for breast implant removal are complications such as capsular contracture, implant rupture and implant malposition, as well

as the patient's desire to change the implant size or shape. The surgeon is advised to use his/her clinical judgment in choosing from current and accepted breast implants removal and replacement surgical techniques to minimize the incidence of adverse reactions and achieve the best results for the patient.

Disposal Method

Treat product that is not returned to the manufacturer as bio hazardous infectious material. Used device can be disposed in suitable disposal unit and subsequently be incinerated by a specialized collection service or in accordance with local regulations.

SPECIFIC INSTRUCTIONS FOR USE

APPLICABLE TO BREAST IMPLANTS CONTAINING A MICROTRANSPONDER

Description

Sterile Silicone Breast Implants Motiva Implant Matrix® with *microtransponder* include a miniaturized, implantable, radio frequency identification device (RFID), which is safely embedded in the breast implant filler material. Scanners to scan and read the information in the *microtransponders* are purchased separately. The *microtransponder* is a passive device that contains an electronic circuit that is activated externally by a low-power electromagnetic field, emitted by a handheld battery powered scanner. The *microtransponder* is used to store an electronic serial number (ESN). The ESN number is used by patient-approved physicians and other health professionals.

Indications

The *microtransponder* is indicated for use as a miniature, long-term implantable device that is inserted into the breast implant. The *microtransponder* in the breast implant provides the patient an electronic serial number that may be used to access a database containing the breast implant information (serial and lot numbers; reference number; volume, size and projection, model, surface type, manufacturing date, etc.).

Contraindications

Breast implants containing a *microtransponder* should not be used in patients known to have allergies or sensitivity to the composition of the *microtransponder* (USP Type III glass).

Precautions

Patients with breast implants containing a *microtransponder* may safely undergo MRI diagnostics in up to 3 Tesla cylindrical systems. See section below: INSTRUCTIONS FOR PATIENTS UNDERGOING MRI for detailed instructions.

Instructions For Patients Undergoing MRI

The patient should be monitored continuously throughout the MRI procedure using visual and audio means (e.g. intercom system). Instruct the patient to alert the MRI system operator of any unusual sensations or problems so that, if necessary, the MRI system operator can immediately terminate the procedure. Provide the patient with a means to alert the MRI system of any unusual sensations or problems.

Do not perform the MRI if the patient is sedated, anesthetized, confused, or otherwise unable to communicate with the MRI system operator.

Patients should be advised to have regular MRIs over their lifetime to screen for silent rupture even if they are not having any apparent problems. As mentioned before, it is advised by the FDA to have the first MRI 3 years postoperatively, then regularly at 2-year intervals.

Sterile Silicone Breast Implants Motiva Implant Matrix® are MRI conditional. The patient implanted with Sterile Silicone Breast Implants Motiva Implant Matrix® can undergo MRI scan under the following conditions:

- Static magnetic field of 1.5-Tesla and 3 -Tesla only.
- Maximum spatial gradient magnetic field of 4.000-gauss/cm (40-T/m) (extrapolated).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode.
- Under the scan defined conditions, the Sterile Silicone Breast Implants Motiva Implant Matrix® with *microtransponder* is expected to produce a maximum temperature rise of 1.5°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the magnetically induced displacement force and magnetically induced torque were tested, and no clinically significant displacement or torque was detected. Sterile Silicone Breast Implants Motiva Implant Matrix® that contains a *microtransponder* creates an imaging void during breast implant MRI (known as artifact effect) that can block visualization of a small area around the *microtransponder*. In non-clinical testing, the image artifact caused by Sterile Silicone Breast Implants Motiva Implant Matrix® extends approximately 15 mm radially from the RFID when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

In selected cases, additional imaging techniques such as ultrasound, tomosynthesis, digital compression mammogram, subtraction contrast mammography and scintimammography are recommended to complement the visualization of the region affected by the artifact and improve the overall diagnosis.

Studies conducted by the manufacturer indicate that the use of a “combined” or “dual” modality, using additional imaging technologies (i.e. MRI with: Ultrasound, Mammography, Tomosynthesis, etc.), may considerably increase the diagnostic accuracy of procedures involving Sterile Silicone Breast Implants Motiva Implant Matrix® with *microtransponder*. The addition of other imaging modalities, using standard practices, allows for the complete radiological survey of the breasts.

Additional Instructions For Use

Additional Directions for Use for breast implants which include a *microtransponder*:

- Verify the *microtransponder* in the implant before opening the sterile barriers, using the corresponding scanner, if available.
- Re-verify the *microtransponder* in the implant after implantation using the corresponding scanner, if available.

Caution

If the breast area is subsequently subjected to physical trauma as a result of accident or injury, the patient shall consult her physician to ensure that the *microtransponder* is properly functioning. If for any reason the *microtransponder* would stop being scannable by the appropriate scanner, this sole situation will not impair the breast implant to continue fulfilling its function appropriately and does not constitute a complication.

SPECIFIC INSTRUCTIONS FOR USE

DEVICE TRACKING

Silicone gel-filled breast implants are subject to device tracking. Compliance with this requirement is mandatory. This means that it is required to report to Establishment Labs, either directly or through a representative, the lot and serial numbers of the device(s) implanted in a patient, the date of the surgery, the ID number and personal contact information, and information related to the surgeon's practice.

Establishment Labs strongly recommends that all patients receiving silicone gel-filled breast implants participate in the Establishment Labs device- tracking program, entering their implant information at www.motivaimplants.com/#/implantRegistration. This will help ensure that Establishment Labs has a record of each patient's contact information so that they can be contacted in the event of a field action or other situations with the implants of which they should be made aware

PROCEDURE RECORDING AND DEVICE ID CARD

Each breast implant is supplied with five Patient Record Labels showing the reference number, lot number, serial number, side (left or right), and volume of the implant. Patient Record Labels are located on the internal product packaging attached to the main label. To complete the Patient ID Card, adhere one Patient Record Label for each implant on the back of each patient's ID Card. Another label should be affixed to the patient's chart. A third label should be attached to the practitioner's records and the fourth label is provided for hospital records when applicable. A spare label is also provided. If a Patient Record Label is unavailable, the appropriate information may be copied by hand from the device label.

PATIENT ID CARD

Every patient must have a record of her surgical procedure in case of future consultations or additional surgeries. Each implant is provided with a Patient ID card, which must be given to the patient for personal reference. Apart from the information stated on the Record Labels that should be affixed to the back of the card, the Patient ID card includes the patient's name, position of the implant (submuscular, subglandular, dual plane, other), date of implantation and name of the treating surgeon.

INFORMATION ON EXPECTED LIFETIME

In practice, it is not possible to predict accurately the actual lifetime of an individual implant. It is well understood that several factors are beyond the control of the manufacture. These factors might have a significant effect on the lifetime of an individual device. These factors include the actual implantation procedure, the anatomy and state of health of the patient and regular activities (for example sporting activities), as well predictable and unpredictable external mechanical influences.

REPORTING AND DEVICE RETURN EFFORTS

In the event of an explantation, the reason for explantation should be reported on the Establishment Labs Complaint Notification Form and the explanted device returned to the local Establishment Labs representative. In case there is no local representative available, report directly at Establishment Labs Coyal Free Zone and Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica; Phone: +506 2434-2400 or <https://motiva.health/support/>.

The explanted device must be decontaminated and properly packaged before returning it according to the Establishment Labs Explanted Implants Return Protocol. In case local health codes do not allow for the implant to be returned, please contact us directly at <https://motiva.health/support/> for specific instructions.

PRODUCT EVALUATION

Establishment Labs requires that any complications resulting from the use of this device be brought to the immediate attention of the company through the Sterile Silicone Breast Implants Motiva Implant Matrix® complaint notification form addressed to Establishment Labs Coyol Free Zone and Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica; Phone: +506 2434-2400 or <https://motiva.health/surgeons/>.

RETURNED GOODS POLICY

Product returns should be handled through the local Establishment Labs representative. In case there is no local representative available, report directly at Establishment Labs Coyol Free Zone and Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica; Phone: +506 2434-2400 or www.motiva.health.

All package seals must be intact for goods to be eligible for return. Returned products may be subject to a restocking charge. For more information, please contact the local Establishment Labs representative.

ESTABLISHMENT LABS ALWAYS CONFIDENT WARRANTY® LIMITED WARRANTY, LIMITATION OF LIABILITY, AND DISCLAIMER OF OTHER WARRANTIES

The complete terms, conditions and limitations of the Establishment Labs Always Confident Warranty® can be reviewed in the website www.motivaimplants.com or can be provided by the local Establishment Labs representative. No warranty or program of Establishment Labs covers any costs, fees or expenses related to any medical treatment and/or the surgical replacement of the implants. Establishment Labs warrants that this product is free of manufacture defects at the time of its shipment. Establishment Labs shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. Establishment Labs's sole responsibility in the event that Establishment Labs determines the product was defective when shipped by Establishment Labs, shall be the replacement of the product. Establishment Labs S.A assumes no further liability. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law, or otherwise, including, but not limited to, any implied warranties of merchantability, suitability for use or performance.

MANUFACTURER

Establishment Labs S.A.

Coyol Free Zone & Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica.

Phone: +506 2434-2400

Establishment Labs S.A.

Coyol Free Zone, 2nd Street, Building B-23, Alajuela, Costa Rica.

Zip code: 20113

www.motiva.health/support

www.motiva.health

www.motivaimplants.com

EUROPEAN REPRESENTATIVE

Emergo Europe:

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands.

EUROPEAN IMPORTER

EDC Motiva BVBA

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Antwerp 2160, Belgium

Non-harmonized symbols in product labeling			
 QTY	Quantity of accessories or medical devices included in the package		Specification of size implant "Base" in centimeters (cm).
	MR conditional, the device can be imaged safely under the tested specifications described in Directions for Use		Specification of size implant "Projection" in centimeters (cm)
	Patient Identification		Health care center or doctor
	Date		Medical device
	Patient information website		Serial Number
	Position of the implant		Unique Device Identifier
	Manufacturer		Catalogue number
	Importer		Country of Manufacture

WARNING: Breast implants have been associated with the development of cancer of the immune system, called Breast Implant Associated Anaplastic Large Cell Lymphoma (BIA-ALCL). More detailed information can be found in the BIA-ALCL section.

Patients with breast implants should have regular clinical follow-ups with their surgeon.

CAUTION:

- Only surgeons with qualified training and certified by their countries' corresponding national medical boards should use this product. This product's use by unqualified practitioners may result in inferior aesthetic outcomes and serious adverse effects.
- Federal (USA) law restricts this device to be used by a certified plastic surgeon.

Directions for Use Sterile Silicone Breast Implants Motiva Implant Matrix® Establishment Labs

INTRODUCTION

This product insert provides an overview of essential information about Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix®, including device description, indications for use, contraindications, warnings, precautions, relevant topics that must be discussed with the patient, adverse events, other reported conditions, returned goods policy, product evaluation, warranty, and medical device reporting.

INTENDED USE

Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® is intended to be used in female patients for the following procedures:

- Breast Reconstruction (primary and revision); to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast anomaly, as well as revision surgery to correct or improve the results of a previous breast reconstruction surgery.

INDICATIONS FOR USE

Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® is indicated in female patients of at least 18 years old for the following conditions:

- To correct congenital or acquired breast deformities or breast asymmetry.
- To correct or improve the result of a previous breast reconstruction surgery.

DEVICE DESCRIPTION AND PERFORMANCE

Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® are invasive, long-term implantable medical devices for breast reconstruction. All implants are composed of a shell, a patch, and highly viscous and highly elastic Establishment Labs silicone gel formula, ProgressiveGel® PLUS or ProgressiveGel® ULTIMA®. The shell is constructed of successive cross-linked layers of silicone elastomer and a low diffusion barrier layer technology that provides these implants their elasticity and integrity. Motiva Implant Matrix® Sterile Silicone Breast Implants are available with and without microtransponder. The microtransponder is placed inside the Motiva Implant Matrix® Sterile Silicone Breast Implants filler material.

Establishment Labs has performed a thorough design verification and validation and non-clinical testing on the Motiva Implant Matrix® Sterile Silicone Breast Implants, including surface characterization, biocompatibility, mechanical, chemical, sterilization, usability, electromagnetic compatibility and electrical safety, and MRI safety studies. Having established the non-clinical safety and performance of the device as it is intended to be used, Establishment Labs continued developing clinical data using the device, supported by input gathered from equivalent devices.

The clinical evidence demonstrates conformity with the relevant safety and performance requirements. The device's performance and safety, as claimed, have been established, and the risks associated with the use of the device are acceptable when weighed against the benefits to the patient.

All raw materials are medical grade, long-term implantable, and biocompatible, as required by international standards.

Basic UDI-DIs of the products covered by this document are:

- Motiva® Round SilkSurface®/SmoothSilk® PLUS: 7445161SilkPlusS4
- Motiva® Round SilkSurface®/SmoothSilk® PLUS with Qid®: 7445161SilkPlusQKU
- Motiva® Round SilkSurface®/SmoothSilk® PLUS with Zen®: 7445161SilkPlusZen8Z
- Motiva Ergonomix® Round SilkSurface®/SmoothSilk®: 7445161ErgoroundR5
- Motiva Ergonomix® Round SilkSurface®/SmoothSilk® with Qid®: 7445161ErgoroundQYN
- Motiva Ergonomix® Round SilkSurface®/SmoothSilk® with Zen®: 7445161ErgoroundZen8Y
- Motiva Round VelvetSurface® PLUS: 7445161VelvetPlus62
- Motiva Round VelvetSurface® PLUS with Qid®: 7445161VelvetPlusQHA
- Motiva Round VelvetSurface® PLUS with Zen®: 7445161VelvetPlusZenEJ
- Motiva Ergonomix® Round VelvetSurface®: 7445161ErgoVelvetM5
- Motiva Ergonomix® Round VelvetSurface® with Qid®: 7445161ErgoVelvetQZ2
- Motiva Ergonomix® Round VelvetSurface® with Zen®: 7445161ErgoVelvetZenPE

The information on the materials to which the patient can be exposed is detailed in the following table.

Motiva Implants® family	Shell			Patch System		Gel	Microtransponder	
	Silicone standard dispersion (% w/w)	Silicone barrier dispersion (% w/w)	Color master-batch (% w/w)	Patch (% w/w)	Dip coat (% w/w)	Silicone filling gel (% w/w)	Ferrite micro-transponder (Qid®) (% w/w)	Non-ferrite microtransponder (Zen®) (% w/w)
Round SmoothSilk®/SilkSurface® Plus	3.09-7.13	0.63-1.46	0.0003-0.0006	0.159-1.572	0.0002-0.0019	89.77-96.11		
Ergonomix® Round SmoothSilk®/SilkSurface®	2.32-5.37	0.58-1.34	0.0002-0.0005	0.159-1.572	0.0002-0.0019	91.66-96.94	0.007-0.066	0.087-0.858
Round VelvetSurface® Plus	2.92-7.43	0.60-1.52	0.0002-0.0006	0.159-1.572	0.0002-0.0019	89.48-96.32		
Ergonomix® Round VelvetSurface®	2.23-5.25	0.56-1.32	0.0002-0.0005	0.159-1.572	0.0002-0.0019	91.85-97.06		

The potential toxicity of the chemicals and metals listed in the following tables have been evaluated with both toxicity testing and risk assessments to assess the exposure levels in comparison to the amount determined to likely be safe. Based on the current results and the risk analysis performed, the leachables/extractables from the shell/patch and gel/ microtransponder of the Ergonomix® Round and Round Plus breasts implant are unlikely to pose a toxicological safety concern.

Quantification of leachable elements digested in concentrated nitric acid by Inductively Coupled Plasma/Mass Spectrometry (ICP/MS).

Element	Concentration (µg/unit)
Barium	0.486
Calcium	27.0
Chromium	1.92
Copper	0.357
Iron	2.13
Magnesium	2.33
Nickel	0.199
Palladium	0.09
Platinum	0.343
Potassium	10.77
Silicon	422.33
Sodium	155.6
Titanium	11.4
Zinc	58.0
ND, not detected	

Extractable Organics Summary for volatile (VOC), semi-volatile (SVOC) and non-volatile (NVOC) compounds in solvents with different polarity indexes: purified water (PW), Hexane (Hex), Ethanol (EtOH), DMC (dichloromethane) and DMSO (dimethyl sulfoxide).


Compound	Concentration (µg/unit)
VOC (HS-GC/MS)	
Trimethyl silanol	415
Benzene	1.03
SVOC (GC/MS)	
Benzoic acid	5,81
Caprolactam	53.6
4-Chlorobenzoic acid	189
4-Chlorobenzoic acid, trimethylsilyl ester	32.7
2,4-Dichlorobenzoic acid	328.9
Decamethyl cyclopentasiloxane (D5)	120.3
Dodecamethyl cyclohexasiloxane (D6)	748.1
Tetradecamethyl cycloheptasiloxane (D7)	513.8
Hexadecamethyl cyclooctasiloxane (D8)	165.5
Octadecamethyl cyclononasiloxane (D9)	391.5
2,2,4,4,6,8-Hexamethyl-6,8- diphenyl- cyclotetrasiloxane	7,506
Eicosamethyl cyclodecasiloxane (D10)	1,053
2,2,4,4,6,6,8,10-Octamethyl-8, 10-diphenyl cyclopentasiloxane	19,485
Cyclic polydimethylsiloxane oligomer (combined values)	113,059
Dimethylsiloxanemethylphenylsiloxane copolymer (combined values)	54,717
Siloxane** (combined values)	11,729.8
1,3,5,7-Tetramethyltetraphenyl cyclotetrasiloxane	19,884
2,4,6,8,10-Pentamethyl-2,4,6,8, 10-pentaphenyl-1,3,5,7,9,2,4,6,8, 10-pentaopentasiloxane isomer (combined values)	50,794
Polymethylphenylsiloxane oligomer (combined values)	23,974
Octamethyl cyclotetrasiloxane (D4)	79.2
2-Ethyl-hexanol	79.3
Tetracosamethyl cyclodecasiloxane (D12)	804
Linear polydimethylsiloxane oligomer (combined values)	377
NVOC (LC/UV)	
Palmitic acid	158.5
Stearic acid	168.2
Erucamide	43.79
Irganox 245	23.7
NVOC (LC/UV-Vis)	
Siloxane** (combined values)	250,375
Di(2-ethylhexyl)phthalate	9,439
Unknown*** (n = 9)	576.87-1,308
GC/MS, gas chromatography/mass spectrometry; LC/UV, liquid chromatography/ultraviolet; LC/UV-Vis, liquid chromatography/ultraviolet visible; NVOC, non-volatile organic compounds; SVOC, semi-volatile organic compounds; VOC, volatile organic compounds.	

Values are reported for the solvent that exhibited higher compound concentration.

The following are the reference ranges of Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix®.

Motiva Implant Matrix® - Round								
Base (cm)	MINI		DEMI		FULL		CORSE	
	P (cm)	V (cc)	P (cm)	V (cc)	P (cm)	V (cc)	P (cm)	V (cc)
8.5	2.2	105	3.1	135	3.5	145	4.0	180
9	2.3	125	3.3	155	3.7	175	4.2	210
9.5	2.4	140	3.4	180	3.9	205	4.5	240
9.75	2.4	150	3.4	190	4.0	220	4.6	260
10	2.5	160	3.5	205	4.1	235	4.8	280
10.25	2.5	170	3.5	215	4.2	255	4.9	300
10.5	2.6	185	3.6	230	4.3	275	5.1	325
10.75	2.6	205	3.7	245	4.4	295	5.2	350
11	2.7	220	3.8	265	4.5	315	5.4	380
11.25	2.7	230	3.8	285	4.6	335	5.5	410
11.5	2.8	245	3.9	300	4.7	355	5.7	440
11.75	2.8	260	3.9	320	4.8	375	5.8	475
12	2.9	275	4.0	340	4.9	400	6.0	510
12.25	2.9	290	4.0	360	5.0	425	6.1	550
12.5	3.0	310	4.1	380	5.1	450	6.3	590
13	3.1	360	4.3	425	5.3	500	6.6	650
13.5	3.2	400	4.4	475	5.5	550	6.9	725
14	3.3	430	4.5	525	5.7	625	7.2	825
14.5	3.4	475	4.6	575	5.9	700	7.5	925
15	3.5	525	4.8	625	6.1	775	7.8	1050

Volumetric Projection

INTENDED CONDITIONS FOR USE

The Motiva Implant Matrix® Sterile Silicone Breast Implants should only be used by licensed physicians/surgeons certified by the corresponding national medical board of their countries, with qualified training in breast implant procedures under sterile conditions, in compliance with good aseptic practices.

INTENDED PATIENTS

The Motiva Implant Matrix® Sterile Silicone Breast Implants are intended for women of at least 18 years of age.

EXPECTED CLINICAL BENEFITS

The following benefits are expected from Motiva Implant Matrix® Sterile Silicone-Gel Breast Implants:

- Reconstruct the breast to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality or
- Revision procedure to correct or improve the results of previous breast reconstruction surgery.

CONTRAINDICATIONS

Establishment Labs Motiva Implants® Sterile Silicone Breast Implants are contraindicated in:

- Women with existing carcinoma of the breast, without mastectomy.
- Women with active infections.

- Women who are currently pregnant or nursing.
- Women with uncontrolled diabetes, which is clinically known to impact wound-healing ability.
- Women who show tissue characteristics clinically incompatible with mammoplasties, such as tissue damage due to compromised vascularity or ulceration.
- Women with any condition – or treatment – determined by the surgeon to constitute an unjustifiable surgical risk (e.g., unstable cardiovascular disease, coagulopathies, chronic pulmonary problems, etc.).

WARNINGS

Care during surgical insertion and subsequent procedures:

- Do not allow sharp objects, such as scalpels or needles, to contact the device during the implantation or other surgical procedures.
- Do not immerse the implant in iodine solution. If iodine solution is used in the pocket, make sure that it is rinsed thoroughly with deionized water so that no residual solution remains in the pocket.
- Do not allow the implant to contact cauterization devices.
- Do not alter the implant or attempt to repair or insert a damaged implant.
- Ensure no excessive force is applied to a very small shell area during the insertion of the device throughout the incision. Instead, apply force over as large an area of the implant as possible during insertion. Excessive force can lead to implant damage by either gel fracture or implant rupture.
- Do not use the periumbilical approach to place the implant.
- Do not place more than one implant per breast pocket.
- Do not treat capsular contracture by closed capsulotomy or forceful external compression, which will likely produce implant damage, rupture, folds, and/or hematoma.
- Do not re-use or re-sterilize any product that has been previously implanted. Breast implants are intended for single use only; therefore, ignoring this indication can cause infection, seroma, capsular contracture, reoperation, unsatisfactory aesthetic results, among other related adverse effects that compromise the patient's health.
- Do not use microwave diathermy in patients with breast implants, as it has been associated with tissue necrosis, skin erosion, and implant extrusion.

PRECAUTIONS

1. Specific Populations

Safety and effectiveness of breast surgery have not been established for the following populations and/or conditions:

- Patients with autoimmune diseases (e.g., lupus, scleroderma).
- Patients whose immune system is compromised (e.g., currently receiving immunosuppressive therapy such as steroids).
- Patients with conditions or medications that may interfere with wound healing ability (e.g., poorly controlled diabetes or corticosteroid therapy) or blood clotting (e.g., concomitant Warfarin therapy).
- Patients with reduced blood supply to breast or overlying tissue.
- Patients who are undergoing radiation therapy.
- Women with ptotic breasts where nipples fall below the inframammary fold, without concurrent mastopexy.
- Previous repeated contour correction failures.
- Patients with a clinical diagnosis of depression or other mental health disorders, including BDD (body dysmorphic disorder) and eating disorders. The patient should be advised to discuss any history of mental health disorders with her surgeon before surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until stabilization of these conditions before undergoing breast implant surgery.
- Other patients with complicated medical histories may be considered to present risk factors that might interfere with the breast implant surgery's safety and effectiveness.

As with any surgical procedure, the patient's medical history should be carefully reviewed to ensure that she is an appropriate candidate for breast implant surgery.

2. Surgical Precautions

Preliminary product examination: Immediately before insertion, examine the device by gently manipulating it while carefully checking for rupture, gel fracture, leakage sites, or particulate contamination.

Surgical technique and implant selection: Several surgical techniques can be used to perform a silicone gel-filled breast implant surgery. Therefore, the surgeon is advised to use his/her clinical judgment to choose the best for the patient, consistent with this product insert.

- The incision should be of appropriate length to accommodate the volume and profile of the implant, reducing the potential for creating excessive stress to the implant when inserting it. Forcing implants through a very small opening may result in local weakening of the breast implant shell, potentially leading to shell damage, gel fracture, and possible implant rupture. After setting realistic aesthetic goals that assure mutual understanding between doctor and patient, the surgeon must choose from current and accepted surgical techniques to minimize the incidence of adverse reactions and achieve the best results.
- Avoid creating wrinkles or folds in the device during the insertion. It is recommended to run a finger around the implant before closing to ensure the implant is flat.
- Procedures such as open capsulotomy, breast pocket revision, hematoma/seroma aspiration, biopsy, and lumpectomy might result in damage to the implant shell, so they must be carefully performed. Care should be taken when repositioning the implant during subsequent procedures to avoid contamination of the implant.
- Periareolar and axillary incision sites may make the insertion more difficult, increasing the risk of damage to the implant. The periareolar incision may considerably reduce the possibility of future breastfeeding.
- The implant size should be consistent with the patient's chest wall dimensions, including base width measurements, characteristics of the tissue, and the implant's projection.
- Textured implants, larger implants, sub-glandular placement, and an insufficient amount of tissue available to cover the implant may cause them to be more palpable.
- Implants of larger sizes may increase the risk of complications such as extrusion, hematoma, infection, palpable implant folds, and visible skin wrinkling.

RELEVANT TOPICS THAT MUST BE DISCUSSED WITH THE PATIENT

Patient Counseling Information

This document and the information for the patient must be carefully reviewed before counseling a patient about Establishment Labs Sterile Silicone Breast Implants Motiva Implant Matrix® and breast surgery. Doctors must thoroughly read and understand the contents of these documents and make sure that any questions or concerns have been resolved before proceeding with the device's use. Breast implant surgery is an elective procedure, and the patient must understand its potential risks and benefits to make an informed decision. For this reason, the patient should be instructed to read the document called "Breast Augmentation and Reconstruction with Motiva Implants® Information for the Patient", which is available on the website <https://ifu.motiva.health/>. The doctor must discuss with the patients the sections on warnings, contraindications, precautions, essential factors to consider, complications², and all other aspects of the document. The physician should inform the patient about the potential complications and that medical management of serious complications may include additional surgery and explanation.

INFORMED CONSENT

It is surgeon's responsibility to document a successful informed decision-making process, providing his/ her "Informed Consent Document", which must be signed by the surgeon, the patient and a witness; this document will become part of the patient's medical record. All possible side effects and complications related to surgery should be advised to patients when considering breast implant surgery with silicone gel breast implants.

As a complement to the informed consent process, Establishment Labs recommends giving each patient the "Breast Augmentation and Reconstruction with Motiva Implants® Information for the Patient" during the surgical consultation; this document is available electronically through the website: <https://ifu.motiva.health/>. Establishment Labs rely on each surgeon's responsibility to ensure that each patient has sufficient time to read and understand the risks, benefits, and recommendations associated with breast implant surgery to make an informed decision.

Summary of Safety and Clinical Performance (SSCP) for the European Community

The SSCP document could be consulted for additional and updated information on the Motiva Implant Matrix® on the company website <https://ifu.motiva.health/> and EUDAMED website <https://ec.europa.eu/tools/eudamed/> using the Basic UDI-DI.

Avoiding Damage During Treatment

Patients should inform other treating physicians of implants' presence to minimize the risk of damage to the implants.

Mental Health and Elective Surgery

All patients seeking an elective procedure such as breast surgery must have realistic expectations that focus on improvement rather than perfection. Ask the patient to openly discuss, before surgery, any history of depression or other mental health disorders.

Breast Examination Techniques

Patients should perform breast self-examinations monthly and be shown how to distinguish the implant from their breast tissue. The patient should not manipulate or squeeze the implant excessively. The patient should be advised that the presence of lumps, persistent pain, swelling, hardening, or change in the implant shape could suggest symptomatic rupture of the implant.

Postoperative Care

The patient should be advised that she will likely feel tired and sore for several days following the operation; her breasts may remain swollen and sensitive to physical contact for a month or longer. She could also experience tightness in the breast area as the skin adjusts to the new breast size. The patient should avoid any strenuous activities for at least a couple of weeks but should be able to return to work within a few days. Breast massage may also be recommended as appropriate. If any issues are present, the patient should report them immediately and possibly have an MRI evaluation to screen for rupture.

Metal Detectors

Motiva Implants® with a microtransponder could be detected by high-sensitivity metal detectors at some airports.

Topical Medications – The patient should consult a physician or a pharmacist before using topical medicines (e.g., steroids) in the breast area.

Smoking – Smoking may interfere with the healing process.

Radiation to the Breast – The literature suggests that radiation therapy may increase the likelihood of capsular contracture, necrosis and implant extrusion.

Insurance Coverage – Patients should check with their insurance company regarding coverage issues before undergoing surgery.

Surgical Setting and Anesthesia – Breast implant surgery is performed in a specialized operating room in a clinic/hospital. General anesthesia is commonly used, and local anesthesia with sedation is also an option. You should be sure to check with your surgeon and with the facility where the surgery will take place to become aware of the tests, pre-surgical examinations, and length of time you need to be without food or your routine medications before the surgical procedure.

POSSIBLE SIDE EFFECTS

If any of the following or other adverse events occur, fill out a complaint notification form, providing all available information on patient, product information, a reason for the complaint, and a summary of the event and send it to <https://motiva.health/surgeons/>.

Because breast implant surgery is more often performed using general anesthesia, it is associated with the same risks as other invasive surgical procedures. After breast implant surgery, patients might experience swelling, hardness, discomfort, itching, bruising, twinges, and pain over the first few weeks. The undesirable side effects identified are detailed below.

Inflammation/Irritation – Breast implants are no different from any foreign material implanted into the human body, and can trigger a host's protective immune reaction. It is a foreign body response demonstrated by redness, swelling, warmth, pain, and or/loss of function.

This foreign body response is universal and ideally removes or otherwise surrounds the "irritant material" with fibrous tissue to prevent unwanted immune consequences.

Gel Fracture – Gel fracture is defined as a fissure or cracks in the implant's gel when excessive intrinsic forces forcibly separate the silicone gel filling. As a result, the implant's shape is irrevocably lost, leading to the need for implant replacement. It can occur with cohesive silicone gel and most frequently due to exposing the implant to excessive compression forces applied to a small shell area during implant insertion. Gel fracture can also occur due to the development of capsular contracture and may result in device distortion.

The incision should be of appropriate length to accommodate the implant's volume and profile with its highly cohesive gel, which will reduce the potential excessive stress that may damage the implant's gel and possibly implant rupture or gel fracture.

Gel fracture can be detected by ultrasound or magnetic resonance imaging (MRI). Most gel fractures are undetectable by physical examination.

Gel Diffusion – Small quantities of silicone may diffuse/bleed through the elastomer envelope of silicone gel-filled implants. The detection of small quantities of silicone in the periprosthetic capsule, axillary lymph nodes, and other distal regions in patients with apparently intact gel-filled implants has been reported in the literature. Some studies on long-term implants have suggested that gel bleed may contribute to capsular contracture and lymphadenopathy development. On the other hand, the evidence demonstrates that gel bleed is a significant contributing factor to capsular contracture. Other local complications are identified even when there are similar or lower complication rates for silicone gel-filled breast implants than for saline-filled breast implants¹.

Redness/Bruising – Bleeding at the time of surgery can cause the skin to change color. It is an expected symptom from the surgery and is likely temporary.

Unsatisfactory Results/Cosmetic Defect – Unsatisfactory results such as stretch marks, visibility and dissatisfaction with the implant volume may occur. Some of these results may cause discomfort. Pre-existing asymmetry may not be entirely correctable by implant surgery. Revision surgery could be indicated to increase patient satisfaction, but this involves additional considerations and risks. Careful preoperative planning and surgical technique can minimize but not always prevent unsatisfactory results.

Pain – Most women undergoing breast implant surgery with a mammary implant will experience post-operative breast and/or chest pain. While this pain usually recedes in most women as they heal after surgery, it can become a chronic problem in other women. Hematoma, migration, infection, implants that are too large or capsular contracture can cause chronic pain. Sudden, severe pain may be associated with implant rupture. The surgeon should instruct the patient to immediately report if there is significant pain or if pain persists.

Skin Rash/Silicone Reaction/Allergy – In general, cutaneous risks with breast implants seem to be low. However, several reports have documented the presence of cutaneous hypersensitivity-like reactions to breast implants, despite the biological compatibility and presumed inertness of their compounds. Topical and systemic medications may relieve the symptoms and lead to a successful resolution. In some cases, implant removal is required for complete symptom relief.

Swelling – Normal postoperative swelling, which peaks about three to five days after surgery, will amplify chest pressure feelings. It is the body's natural response to the trauma of surgery.

Lactation Difficulties – Women who undergo mastectomies and then have breast implant-reconstruction surgeries may not be able to breastfeed on the affected side due to loss of breast tissue and the glands that produce milk.

Heat Sensation – Sterile silicone breast implants Motiva Implant Matrix® with a microtransponder, under MRI scan-defined conditions, can produce a minimal heat sensation.

Hypertrophic Scarring – Scarring is a natural healing process, and it can take time to see improvements. Hypertrophic scars may happen when there is an excessive production of tissue, which forms the scar. Scars may also be caused because the wound takes too long to heal. Biologically, some people tend to be more susceptible to hypertrophic scars due to their genetic makeup².

Bottoming Out – This refers to the inferior displacement of a breast implant that increases the distance between the nipple-areolar complex and the inframammary fold (IMF) after breast-implant surgery. Risk factors reported in the literature are, but not limited to, the lack of quality of pre-existing breast tissue (i.e., thin subcutaneous tissue, defective dermal elements and breast tuberosity), characteristics of the selected breast implant (such as being overly large), IMF dissection, and the type of implant placement during surgery (i.e., submuscular and subglandular planes)³.

The clinical symptoms resulting from a bottomed-out implant include asymmetry, upward-pointing nipples, sagging, palpability, etc. Appropriate surgical planning can mitigate the possible causes of bottoming out. Recommendations include a careful and individual assessment of mammary tissues,

1 Loureno FR, Kikuchi IS, Andreoli Pinto T de J. Silicone Gel Bleed on Breast Implants. Open Biomater J. 2011;3(October):14-17. doi:10.2174/1876502501103010014

2 Berman B, Maderal A, Raphael B. Keloids and Hypertrophic Scars: Pathophysiology, Classification, and Treatment. Dermatol Surg. 2017 Jan;43 Suppl 1:S3-S18. doi: 10.1097/DSS.0000000000000819. PMID: 27347634.

3 Mañero, Ivan M.D.; Montull, Patricia M.D.; Guisantes, Eva M.D. Bottoming Out: A Simple Technique for Correcting Breast Implant Displacement. Plastic and Reconstructive Surgery. December 2009 - Volume 124 - Issue 6 - p 452e-453e

careful implant selection, application of risk-minimizing surgical techniques, and adequate breast support after surgery. Treatments may vary depending on the severity of the complication and range from a simple sub-mammary fixation to the use of additional supporting materials.

Delayed Wound Healing – Some patients may experience a prolonged wound healing time. Smoking causes a decrease in the blood's oxygen levels, directly affecting the surgical wounds' healing process. Delayed wound healing may increase the risk of infection, extrusion and necrosis and may vary depending on the type of surgery or incision.

Capsular Contracture – A capsular contracture pertains to hypertrophic scar tissue investing in a foreign body or surgically implanted device, compromising the aesthetic outcome, resulting in pain, breast deformity, and often necessitating further operations⁴. Detection of breast cancer by mammography may also be challenging. Capsular contracture may be more common following infection, hematoma and seroma, and the chance of it happening may increase over time. Capsular contracture occurs more commonly in patients undergoing revision surgery than in patients undergoing primary implantation surgery. Capsular contracture is the most common complication following implant-based breast surgery and is one of the most common reasons for reoperation.

Capsular contracture is graded into four levels depending on its severity:

- Baker Grade I: The breast is normally soft and looks natural.
- Baker Grade II: The breast is a little firm but looks normal.
- Baker Grade III: The breast is firm and looks abnormal.
- Baker Grade IV: The breast is hard and painful and looks abnormal.

Patients should also be advised that additional surgery might be needed in cases where pain and/or firmness are severe (i.e., Baker Grades III or IV) and that capsular contracture may happen again after additional surgery.

A closed capsulotomy (i.e., external manipulation of the capsule to “pop” the tissue capsule and open it up) is used as a standard procedure for treating capsular contracture; however, most manufacturers, including Establishment Labs, contraindicate it because it can cause implant rupture.

Infection – Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation⁵. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are more challenging to treat than infections in tissue without an implant present. If an infection does not respond to antibiotics, the implant may have to be removed, with replacement occurring only after the infection is resolved.

As with other surgical procedures, Toxic Shock Syndrome (TSS), a life-threatening condition, has been reported in rare instances following breast implant surgery. TSS symptoms occur suddenly and can include high fever (102° F/38.8° C or higher), vomiting, diarrhea, fainting, dizziness, and/or sunburn-like rash. Patients should contact their doctor immediately for diagnosis and treatment if they experience these symptoms⁶.

Seroma – Seroma is an accumulation of fluid that results from tissue inflammation⁷. The etiology of seroma is known in breast surgery and is connected to a hypovascular milieu or trauma following the surgery.

Often, seromas are reabsorbed by the body over several weeks, but needle drainage is sometimes used to remove the fluid⁸. While seromas do not increase breast cancer risk, they sometimes heal with scar tissue or calcifications that can raise a concern about mammograms in the future. Seroma symptoms most often appear a week to 10 days after surgery; the area may feel tender and swollen, with a discrete

lump and redness arising within a day or two. Early seroma formation is defined as periprosthetic fluid accumulation within the first postoperative year, whereas the late form is any moment beyond that time⁹.

In addition to causing pain, a seroma increases the risk of developing an infection in the breast. Depending on the location, it may also increase pressure over the surgical site and sometimes cause wound dehiscence.

Hematoma – A hematoma is a collection of blood within the breast tissue. Hematomas are one of the several complications that may follow breast implant surgery. Symptoms of hematomas generally include swelling, bruising and pain around the incision area¹⁰.

While most hematomas are small and will completely drain on their own and the blood re-absorbs into the body, those patients that are experiencing moderate to severe pain should undergo a follow-up visit. Most hematomas will either resolve on their own or will only require draining. Drains are small surgical tubes that lead out of the breast with a small bulb attached to collect blood and other fluids.

Rupture – Breast implants rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation but is more likely to occur the longer the implant is in place. The following may cause implants to rupture: damage by surgical instruments, implant stress and weakening during implantation, implant age and design, submuscular rather than subglandular location, the occurrence of post-operative hematomas or seromas, folding or wrinkling of the implant shell, excessive force to the chest, trauma, compression during mammographic imaging, and severe capsular contracture¹¹.

Silicone gel-filled implant ruptures are most often silent; this means that most of the time, neither the doctor nor the patient can determine with the physical examination if the implant has a tear or hole in the shell. The integrity of breast implants (and detection of gel fractures and/or silent ruptures) can be evaluated through several techniques. High-resolution ultrasound (HRUS) is widely accepted by healthcare providers and patients for rupture diagnosis. Additionally, the USFDA recommends magnetic resonance imaging (MRI) surveillance with the first MRI performed three years postoperatively and subsequent MRIs performed every two years after that¹². These recommendations may vary from country to country, so please provide the patient with additional guidance based on your country's current care standards. Establishment Labs does not recommend closed capsulotomy for treating capsular contracture because it can cause implant rupture.

Some symptoms may appear, such as lumps surrounding the implant or in the axilla, change or loss of size or shape of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast. These symptoms are not specific to rupture and may also be experienced by patients who have capsular contracture.

Some cases have been reported suggesting that silicone-implant leakage should be considered in the differential diagnosis of eosinophilia¹³.

Extrusion – Breast-implant extrusion, or breast implant exposure, occurs when the breast skin and tissues that are holding the implant fail, causing the implant to protrude through the skin and become exposed. It happens in less than 2% of patients; it is shortly after breast implant surgery or down the road. Breast-implant extrusion can occur for various reasons, such as improper wound healing due to an infection, trauma, too little soft-tissue coverage, oversized implant coupled with too little tissue coverage, or lack of blood supply. Breast-implant extrusion calls for surgery and removal of the implant¹⁴.

4 Headon H, Kasem A, Mokbel K. Capsular Contracture after Breast Augmentation: An Update for Clinical Practice. *Arch Plast Surg.* 2015 Sep;42(5):532-43. doi: 10.5999/aps.2015.42.5.532. Epub 2015 Sep 15. PMID: 26430623; PMCID: PMC4579193.

5 Song JH, Kim YS, Jung BK, Lee DW, Song SY, Roh TS, Lew DH. Salvage of Infected Breast Implants. *Arch Plast Surg.* 2017 Nov;44(6):516-522. doi: 10.5999/aps.2017.01025. Epub 2017 Oct 27. PMID: 29076316; PMCID: PMC5801791.

6 Franchelli S, Pesce M, Baldelli I, Marchese A, Santì P, De Maria A. Analysis of clinical management of infected breast implants and of factors associated to successful breast pocket salvage in infections occurring after breast reconstruction. *Int J Infect Dis.* 2018 Jun;71:67-72. doi: 10.1016/j.ijid.2018.03.019. Epub 2018 Apr 13. PMID: 29660396.

7 Sood A, Kotamarti VS, Therattil PJ, Lee ES. Sclerotherapy for the Management of Seromas: A Systematic Review. *Eplasty.* 2017 Aug 28;17:e25. PMID: 28890747; PMCID: PMC5575675.

8 Gioacchini M, Bottoni M, Grasseti L, Scalise A, Di Benedetto G. A simple, reliable, and inexpensive method for seroma drainage. *Arch Plast Surg.* 2015 May;42(3):361-2. doi: 10.5999/aps.2015.42.3.361. Epub 2015 May 14. PMID: 26015895; PMCID: PMC4439599.

9 Sforza M, Husein R, Atkinson C, Zaccheddu R. Unraveling Factors Influencing Early Seroma Formation in Breast Augmentation Surgery. *Aesthet Surg J.* 2017 Mar 1;37(3):301-307. doi: 10.1093/as/sjw196. PMID: 28207027.

10 FDA. Risks and Complications of Breast Implants. Published 2020. <https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

11 Hillard C, Fowler JD, Barta R, Cunningham B. Silicone breast implant rupture: a review. *Gland Surg.* 2017 Apr;6(2):163-168. doi: 10.21037/gs.2016.09.12. PMID: 28497020; PMCID: PMC5409893.

12 Update on the Safety of Silicone Gel-Filled Breast Implants (2011) - Executive Summary (2011). Accessed March 26, 2021. <https://www.fda.gov/medical-devices/breast-implants/update-safety-silicone-gel-filled-breast-implants-2011-executive-summary>

13 Levenson, Toby; Greenberger, Paul A; Murphy, Robert (1996). Peripheral Blood Eosinophilia, Hyperimmunoglobulinemia A and Fatigue: Possible Complications Following Rupture of Silicone Breast Implants. 77(2), 119–122. doi: 10.1016/S1081-1206(10)63497-7

14 Breast Cancer. Breast Implant Extrusion. Published 2020. <https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after, This is called 'extrusion.'>

Surgical Wound Dehiscence – Surgical wound dehiscence (SWD) is the separation of the margins of a closed surgical incision that has been made in the skin with or without exposure or protrusion of underlying tissue, organs or implants. Separation may occur at single or multiple regions, involve the incision's full length, and affect one or more tissue layers. A dehiscent incision may or may not display clinical signs and symptoms of infection.

Rotation – Anterior/posterior rotation, also called flipping, has been observed more frequently with cohesive gel implants. The flat base of the implant is positioned anteriorly, deforming the breast of the patient. Proper placement and pocket dissection reduce the risk of occurrence¹⁵.

Flipping can be treated with bimanual manipulation in the office and can be done repeatedly in recurrent cases. However, in some cases, revision surgery may be necessary to reduce pocket dimensions.

Literature has reported that the interaction between breast envelopes, the implant's physical characteristics, and pocket dissection could cause malposition. Other theories include the involution of breast tissue.

Regarding implant characteristics, flipping has been associated with the presence or absence of texturing, implant shape/profile, and the gel-filling ratio (i.e., to what degree the implant has been filled). Other factors such as infection, hematoma/seroma, capsular contracture, dissection, surgeon's experience, physical activity, and external manipulation of the implant could potentially contribute to the development of this complication.

Malposition – Malposition of a breast implant is defined as an incorrect placement during surgery or its shifting from its original position. It is also called displacement/lateralization. Malposition has been a frequent event reported due to its multifactorial causes and it can be expected during the lifetime of the device.

The implant's shifting can be produced by trauma, capsular contracture, gravity, or initially wrong placement¹⁶. The surgeon must plan the operation carefully and conduct the surgery with a technique that can minimize but not completely evade the risk of malposition. The risk associated with this event is dissatisfaction with aesthetic outcomes.

Double Capsule – The double capsule refers to the finding of two distinct capsular layers, separated by an inter-capsular space (ICS), around a breast implant. Although rare, double capsules can occur after breast-implant surgery. This condition's etiology is still undefined, but two main hypotheses could explain this complication's development. The first is associated with the periprosthetic capsule's mechanical delamination that creates an intercapsular space (ICS) resulting from fractures after shearing forces are applied between the inner capsule-prosthesis complex and outer capsule. The second hypothesis concerns the development of periprosthetic seromas developing around the capsule. Double capsules may be partial or complete. Clinical signs may vary from asymptomatic to firmness of the implant, discomfort, change in shape or position of the implant, and pain.

Symmastia – This is a relatively rare implant-displacement problem that occurs when the skin and muscle between the breasts over the sternum (breastbone) detaches and the two pockets of tissue that hold the breast implants come together to form one pocket. This allows the implants to come together in the middle, creating the appearance of a "unibob" and sometimes causing discomfort or pain. It's often difficult to correct symmastia and it may require more than one surgery. In most cases, surgery will involve removing the implants and replacing them with new (usually smaller) implants.

Asymmetry – Preoperative asymmetries include areolas in different positions medially or regarding height, different breast shapes (e.g., one round and the other tuberos), or different breast sizes. These asymmetry types should be differentiated from a postoperative aesthetic difference in the two breasts produced by factors described previously, such as a fall of the fold, a high implant, or a rotation of the implant. Asymmetries caused by the implant's unequal fitting or by creating different sub-mammary grooves. They can be prevented using adequate preoperative planning, correct dissection of the pockets, and comparing the two breasts after fitting the implants. It is possible that after breast implant surgery, small deformities in the wall of the thorax or a morphologic mammary disorder may become much more evident. For this reason, the predicted correction of these anomalies should be discussed with the patient before her operation¹⁷.

Rippling/Wrinkling – Rippling is the cutaneous manifestation, visible or palpable, of the implant ripples and edge that are typically most apparent when the patient bends forward. In situations where the implant's soft tissue coverage is insufficient, these harmful effects become more apparent. Risk factors for rippling are related to the breast-tissue quality and low cohesivity of the implanted gel. Adequate coverage over the implant ripples is a mandatory element in preventing rippling or implant edge visibility.

Ptosis – The "waterfall effect" is a descriptive term to indicate sliding ptosis of parenchymal breast tissue over a fixed or encapsulated implant. It occurs more frequently than surgeons anticipate and especially over the long term after augmentation. Certain breast implants are more prone to contribute to this problem as are implants placed in submuscular pockets that ride high, especially in women with anatomical musculoskeletal variance or asymmetry¹⁸.

Calcification/Calcium Deposits – Calcium deposits can form in the scar tissue surrounding the implant and cause pain and firmness. Calcium deposits are visible on mammography. As these deposits must be differentiated from calcium deposits that are a sign of breast cancer, additional surgery may be necessary to remove and examine calcifications. Calcium deposits also occur in women who undergo breast-reduction procedures, in patients who have experienced hematoma(s) and/or seroma(s), and even in the breasts of women who have not undergone any breast surgery. The occurrence of calcium deposits significantly increases with age.

Reoperation (Explantation) – Rupture, unacceptable cosmetic outcomes (i.e., dimpling, wrinkling, and other potentially permanent cosmetic changes of the breast), and other complications may require additional surgeries to the patient's breasts.

Implants are not lifetime devices, and there is a possibility that patients will undergo implant removal(s), with or without replacement, throughout their life. When implants are explanted without replacement, changes to the patient's breasts may be irreversible.

Patients should be advised that the risk of severe capsular contracture doubles for reconstruction patients with implant replacement compared to first-time implantation. There is a risk of an accidental compromise of implant-shell integrity during reoperation, potentially leading to product failure.

Irradiation may cause premature removal because of extrusion, capsular contracture, and recurrent seroma/hematoma.

Visibility/Palpability – Visibility and palpability of the implant can occur because the envelope is thin due to several causes, including excessive volume, the implant's content not being cohesive, or previous surgery and cutaneous aging. If the implant is fitted in a subglandular pocket, a change to a submuscular pocket is needed. The implant volume should be reduced, and it should be ensured that the content is a cohesive gel.

Lymphadenopathy – Silicone-induced lymphadenopathy is a well-known rare complication of implant insertion. It is a disease of the lymph nodes (small, round structures that operate as part of the body's immune system). They become abnormal in size or consistency (most commonly producing swollen or enlarged lymph nodes)¹⁹.

After implant insertion, axillary lymphadenopathy causes are multifactorial, with possible factors including granulomatous reaction, inflammation, and/or malignancy. Literature reports associate lymphadenopathy with intact and ruptured silicone breast implants since microscopic silicone droplets can migrate to body tissues even when the implant surface remains intact. Breast implant rupture and/or leakage through an intact surface can cause fibrosis and granulomatous reactions, which in turn can result in a contracture or regional lymphadenopathy that sometimes mimics malignancy. Various patterns of lymphadenopathy and even extranodal pathology may be observed.

Tissue examination is essential to identify the cause of lymphadenopathy. When in doubt, spectrometry analysis can confirm a diagnosis of silicone-induced lymphadenopathy.

Changes in Nipple and Breast Sensation – Breast surgery can result in an increased/decreased breast and/or nipple sensitivity. Typically, the sensation is lost after complete mastectomy where the nipple itself is removed and can be severely lessened by partial mastectomy. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. While some of these changes can be temporary, they can also be permanent and may affect the patient's sexual response or breastfeeding ability²⁰.

15. Jong, Justin MD*; Gabriel, Allen MD, FACS†; Trekell, Melissa MD†; Lawser, Amy S, MSN, RN§; Heidel, Eric PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T, MD, FACS** Cohesive Round Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321 doi: 10.1097/GOX.0000000000003321

16. Karan Chopra, MD, Arvind U. Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, W. Grant Stevens, MD, Techniques to Repair Implant Malposition after Breast Augmentation: A Review, Aesthetic Surgery Journal, Volume 36, Issue 6, June 2016, Pages 660–671. <https://doi.org/10.1093/asj/sjv261>

17. Breast Asymmetry (2019). Accessed on March 26, 2021. <https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/asymmetry>

18. Frame J. The waterfall effect in breast augmentation. *Gland Surg.* 2017 Apr;6(2):193-202. doi: 10.21037/gls.2016.10.01. PMID: 28497023; PMCID: PMC5409900

19. Lee Y, Song SE, Yoon ES, Bae JW, Jung SP. Extensive silicone lymphadenopathy after breast implant insertion mimicking malignant lymphadenopathy. *Ann Surg Treat Res.* 2017 Dec;93(6):331-335. doi: 10.4174/ast.2017.93.6.331. Epub 2017 Dec 1. PMID: 29250513; PMCID: PMC5729128.

20. Araco A, Araco F, Sorge R, Gravante G. Sensitivity of the nipple-areola complex and areolar pain following aesthetic breast augmentation in a retrospective series of 1200 patients: Periareolar versus submammary incision. *Plast Reconstr Surg.* 2011;128(4):984-989. doi: 10.1097/PRS.0b013e3182268d73

For some patients, additional sensitivity in the nipple area may be noticed days or weeks after breast-implant surgery. This is normal and is due to the stretching of the area near the nerves during surgery. Fortunately, this extra sensitivity will go away as the tissues continue healing. The risk of extra sensitivity is small, and it is typically not permanent.

Necrosis – Necrosis is the formation of dead tissue around the implant. This may prevent wound healing and require surgical correction and/or implant removal. Permanent scar deformity may occur following necrosis. Factors associated with necrosis include infection, steroids in the surgical pocket, smoking, chemotherapy/radiation, and excessive heat or cold therapy.

Interference With Mammography – The importance of mammography should be emphasized. Patients should be instructed to inform their examiners about the presence, type and placement of their implants and request diagnostic mammography rather than screening mammography²¹. Breast implants may complicate the interpretation of mammographic images by obscuring underlying breast tissue and/or compressing overlying tissue.

Even though the presence of breast implants lessens tissue-compression range during mammography, several studies looking at breast cancers in women with implants have found no significant difference in the stage of disease at the time of diagnosis, and prognosis appears to be similar in implanted and nonimplanted patients²². Accredited mammography centers, technicians with experience in imaging patients with breast implants, and displacement techniques are needed to visualize breast tissue in the implanted breast adequately. Anterior breast tissue is best visualized with displacement views and posterior breast tissue with compression views. The decrease in a visible area of 35 % with compression views is improved to 25 % with displacement views.

Interference With Magnetic Resonance Imaging (MRI) – Motiva Implants® with a microtransponder are considered MRI conditional. During the MRI study, the microtransponder can create an MRI artifact immediately around the microtransponder (known as an artifact) that can prevent radiologists from being able to see parts of the implant's footprint and parts of the patient's tissue.

Therefore, there are potential, added MRI risks associated with this artifact including, but not limited to, an inadequate evaluation of the implant shell for the detection of rupture or missing a diagnosis of cancer should it obscure cancer in the artifact area.

The calculated risk of missing a shell rupture due to the artifact is 1 for every 166,000 units of Motiva Implants® with a microtransponder.

The risk of missing breast-cancer detection due to the artifact has been determined to be one high-risk patient with a cancer recurrence for every 596 high-risk patient MRI screening exams performed on patients with Motiva Implants® with a microtransponder. When MRI is used in combination with ultrasound (US) to screen the high-risk patient group, it would take 17,892 MRI and US combination screening exams before a patient with cancer recurrence is likely to be missed (false negative).

Reduction of these risks can be achieved by performing an US in addition to the MRI thereby allowing the radiologist to see the area within the artifact produced in the MRI. As such, the radiologist should be informed of the presence of the microtransponder and that it is embedded near the patch area inside the breast implant. The presence of an MRI imaging void artifact should be anticipated along with its expected size.

OTHER REPORTED CONDITIONS

There have been reports in literature of other conditions in women with silicone gel-filled breast implants. Many of these conditions have been studied to evaluate their potential association with breast implants. No cause and effect relationship has been established between breast implants and the conditions listed below. Furthermore, it is possible that other risks yet unknown may appear, which could be determined to be associated with breast implants in the future.

Neurological Signs and Symptoms – Some women with breast implants have experienced neurological disturbances (e.g., visual symptoms or alterations in sensation, muscle strength, walking, balance, thinking or memory) or diseases (e.g., multiple sclerosis), and they believe those symptoms are related to their implants. However, there is no evidence in the published literature of a causal relationship between breast implants and neurological disease.

21 Sá Dos Reis C, Gremion I, Richi Meystre N. Study of breast implants mammography examinations for identification of suitable image quality criteria. *Insights Imaging*. 2020;Jan 3:11(1):3. doi: 10.1186/s13244-019-0816-5. PMID: 3190684. PMCID: PMC6942083.

22 Helmich LR, Mellemkjaer L, Gunnarsdóttir KA, Tange UB, Krae C, Møller S, McLaughlin JK, Olsen JH. Stage of breast cancer at diagnosis among women with cosmetic breast implants. *Br J Cancer*. 2003;Mar 24:88(6):832-9. doi: 10.1038/sj.bjc.6600819. PMID: 12644818. PMCID: PMC2377078.

Cancer – Breast-cancer reports in medical literature reveal that patients with breast implants are not at a greater risk of developing breast cancer than those without breast implants. Reports in medical literature indicate that breast implants do not significantly delay breast-cancer detection or adversely affect cancer-survival prognosis in implanted women. Some studies even suggest lower rates of breast cancer in women with breast implants.

Breast Mass/Cyst – A breast cyst is a fluid-filled sac that develops within breast tissue. These sacs form when normal fluid producing glands in the breast enlarge or become blocked²³. Depending on the type of breast implant and the implant placement, the mass sometimes occurs in women with thin breast tissue. Breast cysts are usually found on a self-breast exam. When they are small, they often go unnoticed or may instead be seen on a mammogram.

Breast Tissue Atrophy – Breast atrophy may result from aging or the pressure usually exerted by a large breast implant in relation to the patient breast and chest wall size.

Chest Wall Deformity – The breast implant's pressure may cause the breast tissue to thin and shrink (with increased implant visibility and palpability), potentially leading to chest wall deformity. It can occur while implants are still in place or following implant removal without replacement.

Breast Implant Illness (BII) – Over the past several years, the FDA has received new information about systemic symptoms commonly referred to as breast implant illness (BII), which some patients attribute to their implants. Some people with BII also get diagnosed with a specific autoimmune or connective tissue disease²⁴, but many do not.

Researchers are investigating the symptoms to understand their origins better. These symptoms and what causes them are poorly understood. In some cases, removing breast implants without replacement is reported to reverse BII symptoms.

Symptoms can include central nervous system (CNS) disorders (e.g., brain fog, memory loss, tinnitus, vertigo, headaches, blurred vision and migraines); musculoskeletal disorders (e.g., fibromyalgia, muscle pain, hand discoloration, numbness, headaches and migraines); psychological disorders (e.g., anxiety, panic attacks, and feeling of imminent death); immune/inflammatory disorders (e.g., Raynaud syndrome, scleroderma, Hashimoto's thyroiditis, Sjogren's syndrome, autoimmune disease, recurrent infections, rheumatoid arthritis, night sweats, toxic shock, chronic fatigue, dry eye, sudden food intolerance, systemic lupus erythematosus, and multiple sclerosis); as well as anemia and symptoms related to the cardiorespiratory and genitourinary systems.

Connective Tissue Disease (CTD) – Since the early 1990s, nearly a dozen comprehensive systemic reviews have been commissioned by government health ministries in several countries to examine the alleged links between silicone-gel breast implants and systemic diseases²⁵. No hard evidence has been found to support an association between silicone breast implants and CTDs.

Case reports of women with silicone breast implants and CTD include the following symptoms: nervous system alterations (e.g., brain fog, memory loss, blurred vision, migraines, tinnitus); musculoskeletal disease (e.g., muscle/joint pain, fibromyalgia, numbness/tingling in upper and lower limbs, and slow muscle recovery after activity); immune/inflammatory (e.g., Raynaud syndrome, Sjogren's syndrome, Hashimoto's thyroiditis, scleroderma, recurrent/persistent infections, and rheumatoid arthritis); GI/genitourinary (e.g., reduced libido, pancreatitis, urinary tract infection, metallic taste, choking, sudden disease, acid reflux, gastritis, and weight loss/gain); as well as cardiorespiratory and psychological symptoms.

Recent studies still suggest that this association is possible, given that silicone in breast implants can act as a foreign body that can elicit an inflammatory response. Microscopic particles of silicone from the original surgical site have been found far away from it (e.g., in the liver), suggesting a small number of silicone particles detach from the implants and migrate through the lymphatic or circulatory system to other organs. In theory, they could act as adjuvants and start an inflammatory process in joints or activate the immune system and stimulate autoantibodies' production. Nevertheless, no conclusive data is available in this regard²⁶.

23 The NSW Breast Cancer Institute (2006). *Breast cysts - An information guide for patients*

24 FDA. *Risks and Complications of Breast Implants*. Published 2020.

<https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

25 Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. *Ann Intern Med*. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. Epub 2015 Nov 10. PMID: 26550776.

26 Helmich LR, Lipworth L, McLaughlin JK, Friis S. Breast implant rupture and connective tissue disease: A review of the literature. *Plast Reconstr Surg*. 2007;120(7 SUPPL):1162-69. doi: 10.1097/01.prs.0000286684.50274.f2

Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL)^{27, 28}

– BIA-ALCL is a rare type of T-cell lymphoma involving cells of the immune system. The World Health Organization (WHO), in the year 2016, recognized it as a breast implant-associated disease. The exact number of cases remains difficult to determine due to significant limitations in worldwide reporting and a lack of global implant sales data. It has been reported that most data suggests that BIA-ALCL occurs more frequently following implantation of breast implants with textured surfaces rather than those with smooth surfaces. The French National Agency for Medicines and Health Products Safety (ANSM) has requested manufacturers of textured breast implants to perform biocompatibility testing. Establishment Labs has complied with this request. There is a significant body of medical literature relating to breast implants and the risk of developing ALCL. According to the FDA, all the information reviewed as of the date of FDA's March 2017 notice suggests that "women with breast implants have a very low but increased risk of developing ALCL compared to women who do not have breast implants." Most cases of BIA-ALCL are treated by removal of the implant and the capsule surrounding the implant, and some cases have been treated by chemotherapy and radiation.

The following is a consideration from the FDA to the investigators regarding BIA-ALCL: if you have patients with breast implants, you should continue to provide them routine care and support. BIA-ALCL is a very rare condition; when it occurs, it has been identified most frequently in patients undergoing implant-revision operations for late onset, persistent seroma. Because it has generally only been identified in patients with late onset of symptoms, such as pain, lumps, swelling or asymmetry, prophylactic breast-implant removal in patients without symptoms or other abnormality is not recommended. Current recommendations include the steps below.

- Be aware that most confirmed cases of BIA-ALCL have occurred in women with textured breast implants.
- Provide the manufacturers labeling as well as any other educational materials to your patients before surgery and discuss with them the benefits and risks of the different types of implants.

Consider the possibility of BIA-ALCL when you have a patient with late onset, persistent peri-implant seroma. In some cases, patients presented with capsular contracture or masses adjacent to the breast implant. If you have a patient with suspected BIA-ALCL, refer the patient to an appropriate specialist for evaluation. When testing for BIA-ALCL, collect fresh seroma fluid and representative portions of the capsule and send for pathology tests to rule out BIA-ALCL. Diagnostic evaluation should include cytological evaluation of seroma fluid with Wright Giemsa stained smears and cell-block immunohistochemistry testing for cluster of differentiation (CD) and anaplastic lymphoma kinase (ALK) markers. Develop an individualized treatment plan in coordination with the patient's multi-disciplinary care team. Consider current clinical-practice guidelines, such as those from The Plastic Surgery Foundation (PSF) or the National Comprehensive Cancer Network (NCCN) when choosing your treatment approach.

For the latest statistical data on reported cases, refer to [fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alc/](https://www.fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alc/).

Estimated incidence of BIA-ALCL per device surface

In January 2020, the FDA emitted a medical device report for breast implant-associated anaplastic large cell lymphoma (BIA-ALCL). It included 733 U.S. and global reports of BIA-ALCL of which 496 (68 %) were related to textured breast implants and 28 (4 %) to smooth implants. This data is very similar to the data reported in the previous year, which included a total of 673 BIA-ALCL reports of which 385 (67 %) were related to textured breast implants and 26 (5 %) to smooth implants texture, and one had a history of smooth implant and no known textured implant.²⁹

Although the lifetime prevalence of BIA-ALCL was initially estimated to be one in 30,000 women with a textured implant, more recent single-center studies demonstrate incidence rates ranging from 1:355 to 1:559.³⁰⁻³² According to Loch-Wilkinson et al., macro-textured implants carry a significantly higher risk of developing BIA-ALCL. As they describe, it has been shown that textured implants, with their

greater surface area, promote higher levels of bacterial biofilm growth and that this higher bacterial load produces a significant and linear increase in lymphocyte activation³³, increasing the risk of developing BIA-ALCL 14.11 times compared with microtextured implants.³⁴

There are no cases or incidences of BIA-ALCL reported for Motiva Implants®, information is based on literature.

INSTRUCTIONS FOR USE**Single Use**

This product is intended to be used only in one patient for a single procedure. **DO NOT** reuse explanted implants. Reusing a single-use device could expose patients and staff to risks that outweigh the perceived benefits of using such devices. This product is not intended to be reprocessed in any way and/or used again, not even on the same patient. The reuse of single-use devices can affect their safety, performance and effectiveness. Patients can be exposed to unnecessary risks, such as infection and other associated effects, seroma, capsular contracture and possible reoperation. In addition, it is not possible to assure the proper cleaning and decontamination, the presence of residues of cleaning agents, reaction to endotoxins, exposure to another bio-hazards risk, and/or device failure. This practice may also have legal implications that vary according to each jurisdiction.

Product Traceability

The product traceability labels, provided with each device and located within the internal product packaging, provide product specific information and should be attached to the patient's chart for identification purposes. Labels are also available for the Patient Implant Card and the hospital files, if applicable. The surgeon should encourage the patient to participate in the Establishment Labs device-tracking program, entering their implant(s) information at <https://motiva.health/motivalmagine/>.

Traceability helps ensure that Establishment Labs has a record of each patient's contact information to be contacted in the event of a field action or other problems with the implants they should be made aware of.

Sterile Product

The product is sterilized by the manufacturer using the dry heat sterilization method; each sterile silicone breast implant is supplied in a sealed, double sterile barrier primary package. Use standard procedures to maintain sterility during the transfer of the breast implant to the sterile field. Remove the breast implant from its packaging in an aseptic environment using talc-free gloved hands.

The implant's sterility is maintained only if the thermoform packages, including the package seals, are intact.

DO NOT use the product if the thermoform packages or seals have been damaged or have been unintentionally opened before use.

DO NOT re-sterilize the product.

Avoid prolonged exposure to extreme storage conditions. We recommend keeping these devices at room temperature, at atmospheric pressure, in dry conditions, and away from direct sunlight.

DO NOT implant any device that may appear to have particulate contamination, damage or loss of shell integrity. A sterile backup implant must be readily available at the time of surgery.

DO NOT implant any device that may appear to have leaks or scratches.

27 Clemens MW, Nava MB, Rocco N, Miranda RN. Understanding rare adverse sequelae of breast implants: anaplastic large-cell lymphoma, late seromas, and double capsules. *Gland Surg.* 2017 Apr;6(2):169-184. doi:10.21037/gs.2016.11.03. PMID: 28497021; PMCID: PMC5409903.

28 Piubelli ML, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) Is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi:10.1093/asj/sjz19729 E

29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Accessed February 28, 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplantassociated-anaplastic-large-cell-lymphoma>

30 Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plastic and reconstructive surgery.* 2017;139(5):1042-1050. doi:10.1097/PRS.00000000000003282

31 Cordeiro PG, Ghione P, NIA A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery : JPRAS.* 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery.* 2017;140(4):645-654. doi:10.1097/PRS.0000000000003654

34 Matros E, Shamsuddin MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open.* 2021;9(3). doi:10.1097/GOX.00000000000003499

How to Open a Sterile Product Package

DO NOT expose the breast implant to talc, sponges, towels or other contaminants.

1. A non-sterile team member should open the outer package.
2. Remove the inner package and invert it over the sterile field, allowing the sealed inner thermoform package to slide into the field gently.
3. Use the pull-tab to open the lid of the inner thermoform package.
4. Retrieve the breast implant and examine it for any particulate contamination, damage or loss of shell integrity. If satisfactory, return the breast implant to the inner thermoform tray. At this point, you might slightly rinse the implant with a small amount of saline to remove the static and cover the tray with the lid until implantation to prevent contact with airborne and surgical-field particulate contaminants.

If not satisfactory, replace the device with a sterile backup implant.

Note: A sterile backup implant must be readily available at the time of surgery.

Surgical Technique and Implant Selection

Breast implants may be positioned throughout different incisions, in several different tissue planes, described as the pocket location. The selection of this pocket is often a key process in obtaining the desired outcome in individual patients. A pocket that matches the specific conditions of the patient's physical/historical characteristics while maintaining a low-profile risk of complications is crucial. Surgeons should consider the use of published principles for implant selection that minimize the risks of both short and long-term complications.

**Implant Placement Selection**

Submuscular/Subpectoral—Benefits of this technique have been reported to include, but are not limited to, greater soft tissue coverage, lower risk of implant visibility, palpability, capsular contracture and infection, less interference with mammography.

Risks of complications related to implant placement in the submuscular pocket have been reported, but not limited to complex surgical technique, prolonged post-operative recovery, morbidity and higher reported pain; less control of the breast shape, loss of the position, and the inframammary fold resulting in a somewhat unnatural appearance with considerable asymmetry; higher risks of "high-riding" malposition (superior/lateral), "double-bubble" deformity, and constricted lower pole.

Indications for subglandular or subfascial pocket locations may include breast implant primary and revision reconstruction surgeries, patients requiring a pocket exchange from a previous dual plane, submuscular plane, or subfascial plane after revision surgery for a given complication.

Benefits of this pocket location include, but are not limited to, more effortless surgical technique, enhanced control of the shape, the position, inframammary fold of the breast, shorter surgery time, rapid postoperative recovery, less pain and swelling, and lower risks of animation deformities. Risks of complications related to this surgical technique have been reported but are not limited to increased implant visibility, palpability, "bottoming out" malposition, infection, and capsular contracture.

Subfascial—This approach is an alternative implant pocket. When the subfascial dissection is performed, fiber attachments are preserved and could maintain better implant positioning. Additionally, the pectorals fascia favors the creation of a support structure for the upper pole of the implant and avoiding the upper displacement of the implant, which makes the upper pole more natural. Its positive aspects are related to a painless recovery compared to the submuscular position and provision of the implant with extra tissue coverage.

Dual plane—This technique combines varying degrees of soft tissue coverage between the subpectoral placement and the subglandular placement; this enables the surgeon to adjust to the needs of tissue coverage for each patient. The Dual Plane technique is divided into III types. Dual plane I, II, and III.

Indications for dual-plane pocket location may include breast implant primary and revision reconstruction and primary reconstruction (single/two-stage) surgeries.

Benefits of these techniques have been reported, but not limited to, improved soft tissue coverage, less interference with mammography, better lower pole fullness with improved upper and medial pole contour. This technique's reported risks include, but are not limited to, implant visibility, palpability, capsular contracture, some animation deformities.

During the Surgical Procedure

It is advisable to have more than one breast-implant size in the operating room at the time of surgery to allow for flexibility in determining the appropriate size to be used.

- A backup implant should also be available.
- Be advised that excessive force should not be applied to a small area of the shell during device insertion. Instead, apply force over an area of the implant as large as possible.
- The incision should be of appropriate length to accommodate the volume and profile of the implant with a highly cohesive gel, reducing potential, excessive stress to the implant when inserting it. Forcing implants through a very small opening may damage the implant's gel and cause possible ruptures or gel fractures. In the event of gel fracture during implantation, do not insert the implant into the patient's body and replace it with a new one instead.
- Insufficient pocket dissection increases the risk of rupture and incorrect positioning of the implant. A well-defined, dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.
- All devices should be thoroughly assessed for gel fracture, gel bubbles or any other device failures before being inserted into the patient's body during breast-implant surgery. Do not insert an implant with gel fractures at any point; instead, replace the implant with a new one.
- During explantation, surgeons must assess the breast implant's integrity intraoperatively to identify the presence or absence of rupture, gel fracture and gel migration. In the case of a device failure, the implant should be returned to Establishment Labs for revision.
- **DO NOT** use lubricants during placement since they increase the risk of pocket contamination and may affect the tissue-capsule interface.
- **DO NOT** damage the breast implant with sharp surgical instruments, such as needles and scalpels, blunt instruments (such as clamps and forceps), or over-handling and manipulation during introduction into the surgical pocket.
- **DO NOT** use excessive force during breast-implant placement.
- **DO NOT** manipulate the implant for radial expansion, compression or dissection of the pocket.
- **DO NOT** use more than one implant per breast pocket.

Maintaining Hemostasis/Avoiding Fluid Accumulation

The risk of postoperative hematoma and seroma may be reduced by carefully handling the hemostasis during surgery and possibly by the postoperative use of a closed drainage system. Persistent or excessive bleeding must be controlled before implantation. Any postoperative evacuation of hematoma or seroma must be conducted with care to avoid contamination or damage to the breast implant.

SPECIFIC INSTRUCTIONS FOR USE

APPLICABLE TO BREAST IMPLANTS CONTAINING A MICROTRANSPONDER

Description and Use

Motiva Implants® with a microtransponder include a long-term implantable radio frequency identification device (RFID) safely embedded in the breast-implant-filler material. Scanners to read the information in the microtransponders are sold separately. The microtransponder is a passive RFID that uses radio waves to provide an electronic serial number (ESN) that assures full traceability to implant-specific data.

The microtransponder in the breast implant provides the patient an ESN used to access a database containing the breast-implant information (i.e., serial and lot numbers; the reference number; volume, size, and projection, model, surface type, manufacturing date, etc.).

Precautions

Patients with breast implants containing a microtransponder may safely undergo MRI diagnostics in up to three Tesla cylindrical systems. See the next section called Instructions for Patients Undergoing MRI for detailed instructions.

Instructions for Patients Undergoing MRI

Patients should be advised to have regular MRIs over their lifetime to screen for silent rupture even if they are not having any apparent problems. As mentioned before, it is advised by the FDA to have the first MRI three years post-operatively, then regularly at two-year intervals.

The patient should be monitored continuously throughout the MRI procedure using visual and audio means (e.g., through an intercom system). Instruct the patient to alert the MRI system operator of any unusual sensations or problems so that, if necessary, the MRI system operator can immediately terminate the procedure. Provide the patient with a means to alert the MRI system operator of any unusual sensations or problems.

Do not perform the MRI if the patient is sedated, anesthetized, confused, or otherwise unable to communicate with the MRI system operator.

Motiva Implants® are labeled MRI conditional. The patient implanted with Motiva Implant Matrix® sterile silicone breast implants can undergo an MRI scan under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla only
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m) (extrapolated)
- Maximum MR system reported whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the normal Operating mode
- Under the scan-defined conditions, the Motiva Implant® Matrix® sterile silicone breast implants with a microtransponder are expected to produce a maximum temperature rise of 1.5°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the magnetically induced displacement and magnetically induced torques were tested, resulting in no clinically significant displacement or torque detected. Motiva Implants® that contain a microtransponder create an imaging void during breast-implant MRI (known as artifact effect) that can block visualization of a small area around the microtransponder. In non-clinical testing, the image artifact caused by Motiva Implants® extends approximately 15 mm radially from the RFID when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

In selected cases, additional imaging techniques such as ultrasound, tomosynthesis, digital compression mammogram, subtraction contrast mammography and scintimammography are recommended to complement the region's visualization affected by the artifact and improve the overall diagnosis.

Studies conducted by the manufacturer indicate that the use of a "combined" or "dual" modality, using additional imaging technologies (i.e., MRI with: ultrasound, mammography, tomosynthesis, etc.), may considerably increase the diagnostic accuracy of procedures involving Motiva Implant® Matrix® sterile silicone breast implants with a microtransponder. The addition of other imaging modalities, using standard practices, allows for the complete radiological survey of the breasts.

ADDITIONAL INSTRUCTIONS FOR USE

Here are additional directions for use of breast implants that include a microtransponder:

- Verify the microtransponder in the implant before opening the sterile barriers using the corresponding scanner, if available.

Caution

If the breast area is subsequently subjected to physical trauma as a result of accident or injury, the patient shall consult her physician to ensure that the microtransponder is appropriately functioning. If for any reason the microtransponder were to stop being scannable by the appropriate scanner, this sole situation will not impair the breast implant to continue fulfilling its function appropriately and does not constitute a complication.

Device Tracking

Silicone gel-filled breast implants are subject to device tracking. Compliance with this requirement is mandatory; therefore, it is required to report to Establishment Labs, either directly or through a representative, the lot and serial numbers of the device(s) implanted in a patient, the date of the surgery, the ID number and personal contact information, and the information related to the surgeon's practice.

Establishment Labs strongly recommends that all patients receiving silicone gel-filled breast implants participate in the Establishment Labs device-tracking program, entering their implant information at motiva.health/motivalmagine/. This will help ensure that Establishment Labs has a record of each patient's contact information to be contacted in the event of a field action or other situations with the implants of which they should be made aware.

Instructions and Precautions for Removal

Among the most common reasons for breast-implant removal are complications, such as capsular contracture, implant rupture and implant malposition, and the patient's desire to change the implant size or shape. The surgeon is advised to use his/her clinical judgment in choosing from current and accepted breast-implant removal and replacement surgical techniques to minimize the incidence of adverse reactions and achieve the best results for the patient.

Storage and Handling Conditions

There are no specific identified conditions necessary for the transportation of devices; however, it is recommended that packaged units should be stored in a designated clean area within the hospital/clinic, stored flat under room conditions.

Expiration

The package labeling contains the expiration date of the device. **WARNING:** Establishment Labs does not guarantee sterility if the device is used after the expiration date.

Disposal Method

Treat any product that is not returned to the manufacturer as bio-hazardous and infectious material. Used devices can be disposed of in a suitable disposal unit and subsequently be incinerated by a specialized collection service or following local regulations.

Mishandling of explanted devices can lead to unnecessary risks, such as infection or microbial hazards.

Information on the Expected Lifetime

Motiva Implants® are not lifetime devices. Based on the chemical characteristics of the materials used in Motiva® devices, the accelerated aging testing for five-years shelf-life, PMS information, and a vast literature review from clinical data obtained from equivalent devices, a survival rate at 10 years of 80 % has been established as acceptance criteria for the Motiva Implants® lifespan³⁵.

It is not possible to predict accurately the actual lifetime of an individual implant. It is well-understood that several factors are beyond the control of the manufacturer having a significant effect on the lifespan of an individual device. These factors include the actual implantation procedure, the anatomy, the state of the patient's health, and behavior(s) and regular activities of the patient (for example, sporting activities), as well as predictable and unpredictable external mechanical influences.

Procedure for Recording and the Device Implant ID Card

Each breast implant is supplied with 10 patient record labels showing the reference number, serial number and implant volume besides other regulatory information. Patient record labels are located on the internal product packaging attached to the main label. To complete the patient implant ID card, adhere one patient record label for each implant to the back of each patient's ID card. Another label should be affixed to the patient's chart. A third label should be attached to the practitioner's records. The fourth label is provided for hospital records; when applicable, additional labels could be used as required to assure device traceability. If a patient record label is unavailable, the appropriate information may be copied by hand from the device label.

Patient Implant ID Card

Every patient must have a record of the surgical procedure in case of future consultations or additional surgeries. Each implant is provided with a patient implant card, which must be given to the patient for personal reference. Record labels should be affixed to the patient implant card. The following information is included in the patient implant card: device type, patient's name, the position of the implant (i.e., submuscular, subglandular, dual plane or subfacial plane), date of implantation, name of the treating surgeon, manufacturer information, UDI, device reference, serial number, and implant volume.

PRODUCT EVALUATION

Establishment Labs requires that any complications resulting from the use of this device be brought to the immediate attention of the company through the complaint notification form addressed to Establishment Labs S.A., Coyol Free Zone and Business Park Building, 4th Street, Building B-15, Alajuela, Costa Rica. For additional assistance, call +506 2434-2400 or go to <https://motiva.health/surgeons/>.

REPORTING AND DEVICE RETURN EFFORTS

In the event of an explantation, the reason for explantation should be reported on the complaint notification form and the explanted device returned to the local Establishment Labs representative. In case there is no local representative available, report directly at Establishment Labs S.A., Coyol Free Zone and Business Park Building, 4th Street, Building B-15, Alajuela, Costa Rica. For additional assistance, call+506 2434-2400 or go to motiva.health/patients-support/.

According to the Establishment Labs explanted implants return protocol, the explanted device must be decontaminated and properly packaged before returning it. In case local health codes do not allow for the implant to be returned, please contact us directly at motiva.health/patients-support/ for specific instructions.

Any serious incident related to Motiva® implant devices should be reported to Establishment Labs and the corresponding competent authority in which the patient is established, if applicable according to local regulations.

RETURNED GOODS POLICY

Product returns should be handled through the local Establishment Labs representative. In case there is no local representative available, report directly at Establishment Labs S.A., Coyol Free Zone and Business Park Building, 4th Street, Building B-15, Alajuela, Costa Rica. For additional assistance, call +506 2434-2400 or go to motiva.health.

All package seals must be intact for goods to be eligible for return. Returned products may be subject to a restocking charge. For more information, please contact the local Establishment Labs representative.

PRODUCT ORDERING

To order directly or for product information, please contact your local Establishment Labs representative. In case there is no local representative available, contact Establishment Labs S.A., Coyol Free Zone and Business Park Building, 4th Street, Building B-15, Alajuela, Costa Rica. For additional assistance, call +506 2434-2400 or email customerservice@establishmentlabs.com. For detailed indications, contraindications, warnings and precautions associated with the use of Motiva Implant Matrix® sterile silicone breast implants, please refer to the product insert provided with each product.

ESTABLISHMENT LABS ALWAYS CONFIDENT WARRANTY® LIMITED WARRANTY, LIMITATION OF LIABILITY, AND DISCLAIMER OF OTHER WARRANTIES

The complete terms, conditions and limitations of the Establishment Labs Always Confident Warranty® can be reviewed on the website motiva.health or can be provided by the local Establishment Labs representative. No warranty or program of Establishment Labs covers any costs, fees or expenses related to any medical treatment and/or the surgical replacement of the implants. Establishment Labs shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. In the event that Establishment Labs determines the product was defective when shipped by Establishment Labs, it shall be the sole responsibility of Establishment Labs to replace the product. Establishment Labs S.A. assumes no further liability. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law, or otherwise, including, but not limited to, any implied warranties of merchantability, suitable for use or performance.

LEGAL MANUFACTURER

Establishment Labs S.A.:

Coyol Free Zone & Business Park Building 4th Street, Building B-15,

Alajuela, Costa Rica Zip code: 20113

Phone: +506 2434-2400

www.establishmentlabs.com

motiva.health/patients-support/

<https://ifu.motiva.health/>

MANUFACTURING SITES

Coyol Free Zone & Business Park Building 4th Street, Building B-15,

Alajuela, Costa Rica.

Zip code: 20113

Coyol Free Zone & Business Park Building 0 Street, Building B-25,

Alajuela, Costa Rica

Zip code: 20113

Establishment Labs S.A.

Coyol Free Zone, 2nd Street, Building B-23,

Alajuela, Costa Rica.

Zip code: 20113

EUROPEAN REPRESENTATIVE

Emergo Europe:

Westervootsedijk 60, 6827 AT Arnhem, The Netherlands


EUROPEAN IMPORTER

EDC Motiva BVBA

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Antwerp 2160, Belgium

Symbols used in product labeling (including non-harmonized symbols)

	Quantity of accessories or medical devices included in the package		Specification of size implant "Base" in centimeters (cm)		Keep away from sunlight
	MR conditional, the device can be imaged safely under the tested specifications described in Directions for Use		Specification of size implant "Projection" in centimeters (cm)		Keep dry
	Patient Identification		Health care center or doctor		Do not re-use
	Date		Medical device		Consult instructions for use
	Position of the implant		Serial Number		Caution
	Manufacturer		Unique Device Identifier		Authorized Representative in the European Union
	Importer		Catalogue number		Use-by date
	Double sterile barrier system		Country of Manufacture		Do not re-sterilize
	Sterilized using steam or dry heat		Date of Manufacture		Fragile, handle with care
	Imaging studies conditional		Do not use if package is damaged and consult instructions for use		Refer to instruction manual/booklet
	Specification of size implant "Height" in centimeters (cm) (Note: Not applicable to Motiva® Ergonomix® Round and Motiva® Round PLUS.)		Implant Volume	A single-digit that could be a letter or a number provided in the product label and in the traceability label	On-line registration validation digit

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Update on the Safety of Silicone Gel-Filled Breast Implants (2011) – Executive Summary (2011). 12 Accessed March 26, 2021. https://www.fda.gov/medical-devices/breast-implants/update-safety-silicone-gel-filled-breast-implants-2011-executive-summary

Levenson, Toby; Greenberger, Paul A; Murphy, Robert (1996). Peripheral Blood Eosinophilia. 13 Hyperimmunoglobulinemia A and Fatigue: Possible Complications Following Rupture of Silicone Breast Implants. 77(2), 119–122. doi: 10.1016/S1081-1206(10)63497-7

Breast Cancer. Breast Implant Extrusion. Published 2020. 14 https://www.breastcancer.org/treatment/surgery/breast-reconstruction/corrective-reconstruction/implant-extrusion. In relationly rare cases after extrusion.”

Jong, Justin MD*; Gabriel, Allen MD; FACS; Trekkell, Melissa MD†; Lawser, Amy S. MSN, RN§; Heidel, Eric 15 PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T. MD, FACS* Cohesive Rupture of Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321. doi: 10.1097/GOX.0000000000000321

Karan Chopra, MD, Arvind U. Gowda, MD, Michele Eagan, MD, W. Grant Stevens, MD, 16 Techniques to Repair Implant Malposition after Breast Augmentation: A Review, Aesthetic Surgery Journal, Volume 36, Issue 6, June 2016, Pages 606–671, https://doi.org/10.1093/asj/sjv261

Breast Asymmetry (2019). Accessed on March 26, 2021. https://www.breastcancer.org/treatment/surgery/ 17

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Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma 28 (ALCL) is Distinct from Systemic ALCL. *ALK Negative in a Patient with Gluteal Implants. Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E29

Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. 29 Accessed February 28, 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplantassociated-anaplastic-large-cell-lymphoma>

Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell 30 Lymphoma. *Plastic and reconstructive surgery.* 2017;139(5):1042-1050. doi:10.1097/PPS.0000000000003282

Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma 31 (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery* . JPRAS. 2020;73(5):841-846. doi: 10.1016/j.bjps.2019.11.064

Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: 32 Determining an Accurate Risk. *Annals of Surgery.* 2020;272(3):403-409. doi:10.1097/SLA.000000000000179

Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in 33 Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery.* 2017;140(4):645-654. doi:10.1097/PPS.0000000000003654

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25. Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. Ann Intern Med. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. Epub 2015. Nov 10. PMID: 26550776.

26. Hölmiş LR, Lipworth L, McLaughlin JK, Friis S. Breast implant rupture and connective tissue disease: A review of the literature. Plast Reconstr Surg. 2007;120(7 SUPPL. 1):62-69. doi: 10.1097/01.prs.0000286664.50274.f2

22. Hölmiş LR, Mellemkjaer L, Gunnarsdóttir KA, Tange UB, Krag C, Møller S, McLaughlin JK, Olsen JH. Stage of breast cancer at diagnosis among women with cosmetic breast implants. Br J Cancer. 2003 Mar 24;88(6):832-8. doi: 10.1038/sj.bjc.6600819. PMID: 12644818. PMID: PM2377078.

23. The NSW Breast Cancer Institute (2006). Breast cysts – An information guide for patients

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За последните статистически данни за докладваните случаи вижте [fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alcl](https://www.fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alcl).

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28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL/ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Отворено на 28 февруари 2022. <<https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplant-associated-anaplastic-large-cell-lymphoma>

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery: JPRAAS.* 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

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34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Silicone Implants Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open.* 2021;9(3). doi:10.1097/GOX.00000000000003499.

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6 Franchelli S, Pesce M, Baldelli I, Marchese A, Santi P, De Maria A. Analysis of clinical management of infected breast implants and of factors associated to successful breast pocket salvage in infections occurring after breast reconstruction. *Int J Infect Dis*. 2018 Jun;71:67-72. doi: 10.1016/j.ijid.2018.03.019. Epub 2018 Apr 13. PMID: 29660396.

7 Sood A, Kotamarti VS, Therattil PJ, Lee ES. Sclerotherapy for the Management of Seromas: A Systematic Review. *Eplasty*. 2017 Aug 28;17:e25. PMID: 28890747; PMCID: PMC5575675.

8 Gioacchini M, Bottoni M, Grasseti L, Scalise A, Di Benedetto G. A simple, reliable, and inexpensive method for seroma drainage. *Arch Plast Surg*. 2015 May;42(3):361-2. doi: 10.5999/aps.2015.42.3.361. Epub 2015 May 14. PMID: 26015895; PMCID: PMC4439599.

9 Sforza M, Husein R, Atkinson C, Zaccachedu R. Unraveling Factors Influencing Early Seroma Formation in Breast Augmentation Surgery. *Aesthet Surg J*. 2017 Mar 1;37(3):301-307. doi: 10.1093/asj/sjw196. PMID: 28207027.

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14 Breast Cancer. Breast Implant Extrusion. Published 2020. [https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after. This is called "extrusion"](https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after. This is called).

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15 Jong, Justin MD*; Gabriel, Allen MD, FACS†; Trekkel, Melissa MD‡; Lawser, Amy S, MSN, RN§; Heidel, Eric PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T, MD, FACS** Cohesive Round Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open; December 2020 - Volume 8 - Issue 12 - p e3321 doi: 10.1097/GOX.00000000000003321

16. Karan Chopra, MD, Arvind U. Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, W. Grant Stevens, MD, Techniques to Repair Implant Malposition after Breast Augmentation: A Review, Aesthetic Surgery Journal, Volume 36, Issue 6, June 2016, Pages 660–671, <https://doi.org/10.1093/asj/sjv261>

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23 The NSW Breast Cancer Institute. (2006). *Breast cysts - An information guide for patients*

24 FDA. *Risks and Complications of Breast Implants*. Published 2020. <https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

25 Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants. *American Review*. Ann Intern Med. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. Epub 2015 Nov 10. PMID: 26550776.

26 Hølmich LR, Lipworth L, McLaughlin JK, Fris S. Breast implant rupture and connective tissue disease: A review of the literature. *Plast Reconstr Surg*. 2007;120(7 SUPPL 1):62-69. doi: 10.1097/01.prs.0000286664.50274.12

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28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Accessed February 28, 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplant-associated-anaplastic-large-cell-lymphoma>

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery: JPRAS*. 2020;73(5):841-846. doi: 10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery*. 2020;272(3):403-409. doi: 10.1097/SLA.0000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery*. 2017;140(4):645-654. doi: 10.1097/PRS.00000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgery Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi: 10.1097/GOX.00000000000003499

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23 The NSW Breast Cancer Institute (2006). Breast cysts - An information guide for patients

24 FDA. Risks and Complications of Breast Implants. Publiceret 2020.

https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants

25 Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. *Ann Intern Med*. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. E-publikation 10. nov. 2015. PMID: 26550776.

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28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALCL Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Tilgæet 28. februar 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplant-associated-anaplastic-large-cell-lymphoma>

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery : JPRA.* 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and Reconstructive Surgery.* 2017;140(4):645-654. doi: 10.1097/PRS.00000000000003654

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23 The NSW Breast Cancer Institute (2006). Breast cysts – An information guide for patients

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28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441–NP442. doi: 10.1093/asj/sjz13729 E29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma FDA. Zugriff am 28. Februar 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplantassociated-anaplastic-large-cell-lymphoma>

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31 Cordeiro PG, Ghilone P, NIA et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery* : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery*. 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery*. 2017;140(4):645-654. doi:10.1097/PRS.0000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiological Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi:10.1097/GOX.00000000000003499.

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³¹ Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery* : JPRAS. 2020; 73(5):841-846. doi:10.1016/j.bjps.2019.11.064

³² Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery*. 2020; 272(3):403-409. doi:10.1097/SLA.00000000000004179

³³ Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery*. 2017; 140(4):645-654. doi:10.1097/PRS.00000000000003654

³⁴ Matros E, Shamsuddin MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgery Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021; 9(3). doi:10.1097/GOX.00000000000003499

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28. Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19279 E

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31 Cordeiro PG, Ghione P, Ni A, *et al.* Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of Plastic, Reconstructive & Aesthetic Surgery* : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, *et al.* Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of Surgery*. 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, *et al.* Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and Reconstructive Surgery*. 2017;140(4):645-654. doi:10.1097/PRS.0000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, *et al.* Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi:10.1097/GOX.00000000000003499

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14 Breast Cancer. Breast Implant Extrusion. Published 2020. <https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after This is called 'extrusion.'>

15 Jong, Justin MD; Gabriel, Allen MD, FACS; Trekkell, Melissa MD; Lawser, Amy S, MSN, RN; Heide, Eric PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T, MD, FACS** Cohesive Round Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321 doi: 10.1097/GOX.0000000000003321

16 Khan Chagra, MD, Arvind U, Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, W, Grant Stevens, MD. Techniques to Repair Implant Malposition after Breast Augmentation: A Review. Aesthetic Surgery Journal. Volume 36, Issue 6, June 2016, Pages 660–671. <https://doi.org/10.1093/asj/sjv261>

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²⁸ Plubeli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

Najnovije statističke podatke o prijavljenim slučajevima potražite na poveznici [fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alc](https://www.fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alc).

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²⁹ Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Accessed February 28, 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplant-associated-anaplastic-large-cell-lymphoma>

³⁰ Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plastic and reconstructive surgery.* 2017;139(5):1042-1050. doi: 10.1097/PRS.0000000000003282

³¹ Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery : JPRAS.* 2020;73(5):841-846. doi: 10.1016/j.bjps.2019.11.064

³² Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. doi: 10.1097/SLA.0000000000004179

³³ Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery.* 2017;140(4):645-654. doi: 10.1097/PRS.0000000000003654

³⁴ Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open.* 2021;9(3). doi: 10.1097/GOX.00000000000003499

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14 Breast Cancer. Breast Implant Extrusion. Publikálási dátuma:2020. <https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion> In relatively rare cases after this is called "extrusion."

15 Jong, Justin MD*; Gabriel, Allen MD, FACS†; Trekkell, Melissa MD†; Lawser, Amy S, MSN, RNS, Heidiel, Eric PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T, MD, FACS** Cohesive Round Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321 doi: 10.1097/GOX.00000000000003321

16 Karan Chopra, MD, Arvind U. Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, W Grant Stevens, MD, Techniques to Repair Implant Malposition after Breast Augmentation: A Review. *Aesthetic Surgery Journal*. Volume 36, Issue 6, June 2016. Pages 660–671. <https://doi.org/10.1093/asj/sjv261>

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²³ The NSW Breast Cancer Institute (2006). Breast cysts - An information guide for patients

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²⁵ Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. Ann Intern Med. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. Epub 2015 Nov 10. PMID: 26550776.

²⁶ Hölmlch LR, Lipworth L, McLughlin JK, Friis S. Breast implant rupture and connective tissue disease: A review of the literature. Plast Reconstr Surg. 2007;120(7 SUPPL. 1):62-69. doi: 10.1097/01.prs.0000286664.50274.f2

²⁷ Clemens MW, Nava MB, Rocco N, Miranda RN. Understanding rare adverse sequelae of breast implants: anaplastic large-cell lymphoma, late seromas, and double capsules. Gland Surg. 2017 Apr;6(2):169-184. doi: 10.21037/gls.2016.11.03. PMID: 28497021; PMCID: PMC5409903.

²⁸ Plubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) Is Distinct from Systemic ALCL/ALK Negative in a Patient with Gluteal Implants. Aesthetic Surg J. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

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³⁰ Doren EL, Miranda RN, Seiber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. Plastic and reconstructive surgery. 2017;139(5):1042-1050. doi: 10.1097/PRS.00000000000003282

³¹ Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. Journal of plastic, reconstructive & aesthetic surgery. JPRAS. 2020;73(5):841-846. doi: 10.1016/j.jbps.2019.11.064

³² Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. Annals of surgery. 2020;272(3):403-409. doi: 10.1097/SLA.0000000000004179

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10 FDA. Risks and Complications of Breast Implants. Pubblicato 2020. <https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

11 Hillard C, Fowler JD, Barta R, Cunningham B. Silicone breast implant rupture: a review. *Gland Surg.* 2017 Apr;6(2):163-168. doi: 10.21037/gls.2016.09.12. PMID: 28497020; PMCID: PMC5409893.

12 Update on the Safety of Silicone Gel-Filled Breast Implants (2011) - Executive Summary (2011). Data di consultazione: 26 marzo 2021. <https://www.fda.gov/medical-devices/breast-implants/update-safety-silicone-gel-filled-breast-implants-2011-executive-summary>

13 Levenson, Toby; Greenberger, Paul A; Murphy, Robert (1996). Peripheral Blood Eosinophilia, Hyperimmunoglobulinemia A and Fatigue: Possible Complications Following Rupture of Silicone Breast Implants. *77(2)*, 119–122. doi: 10.1016/S1081-1206(10)63497-7

14 Breast Cancer. Breast Implant Extrusion. Pubblicato 2020. [https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after,This is called "extrusion."](https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after,This is called)

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24 FDA. *Risks and Complications of Breast Implants*. Pubblicato 2020. <https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

25 Balk EM, Earley A, Avendano EA, Raman G. *Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review*. *Ann Intern Med*. 2016 Feb 2; 164(3):164-75. doi: 10.7326/M15-1169. Epub 2015. Nov 10. PMID: 26550776.

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28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. *Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALCL Negative in a Patient with Gluteal Implants*. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

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Per i dati statistici più recenti sui casi segnalati, consultare il documento [fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alc](https://www.fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alc).

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. Journal of plastic, reconstructive & aesthetic surgery : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. Annals of surgery. 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. Plastic and reconstructive surgery. 2017;140(4):645-654. doi:10.1097/PRS.00000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. Plastic and Reconstructive Surgery Global Open. 2021;9(3). doi:10.1097/GOX.00000000000003499.

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. Journal of plastic, reconstructive & aesthetic surgery. JPRAS. 2020;73(5):841-846. doi:10.1016/j.BJPS.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. Annals of surgery. 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

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34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgery Database: Epidemiologic Implications for BIA-ALCL. Plastic and Reconstructive Surgery Global Open. 2021;9(3). doi:10.1097/GOX.00000000000003499

27 Clemens MW, Nava MB, Rocco N, Miranda RN. Understanding rare adverse sequelae of breast implants: anaplastic large-cell lymphoma, late seromas, and double capsules. Gland Surg. 2017 Apr;6(2):169-184. doi:10.21037/gs.2016.11.03. PMID: 28497201; PMCID: PMC5409903.

28 Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) Is Distinct from Systemic ALCL/ALK Negative in a Patient with Gluteal Implants. Aesthetic Surg J. 2019;39(10):NP441-NP442. doi:10.1093/asj/sjz19729 E

29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Piekļūts 2022. gada 28. februārī. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breastimplant-associated-anaplastic-large-cell-lymphoma>

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14 Breast Cancer. Breast Implant Extrusion. Gepubliceerd 2020. [https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after. This is called "extrusion."](https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In%20relatively%20rare%20cases%20after%20this%20is%20called%20%22extrusion%22.)

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23 The NSW Breast Cancer Institute. (2006). Breast cysts - An information guide for patients

24 FDA Risks and Complications of Breast Implants. Gepubliceerd 2020. <https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

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25 Balk EM, Earley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. *Ann Intern Med.* 2 feb 2016; 164(3):164–175. doi: 10.7326/M15-1169. E-pub 10 nov 2015. PMID: 26550776.

26 Hölmlch LR, Lipworth L, McLaughlin JK, Friis S. Breast implant rupture and connective tissue disease: A review of the literature. *Plast Reconstr Surg.* 2007;120(7 SUPPL. 1):62–69. doi: 10.1097/01.prs.0000286664.50274.f2

27 Clemens MW, Nava MB, Rocco N, Miranda RN. Understanding rare adverse sequelae of breast implants: anaplastic large-cell lymphoma, late seromas, and double capsules. *Gland Surg.* apr 2017;6(2):169–184. doi: 10.21037/gs.2016.11.03. PMID: 28497021; PMCID: PMC5409903.

28 Piuelli ML, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL. ALCL Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441–NP442. doi: 10.1093/asj/sjz19729 E

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33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. Plastic and reconstructive surgery. 2017;140(4):645-654. doi:10.1097/PRS.0000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. Plastic and Reconstructive Surgery Global Open. 2021;9(3). doi:10.1097/GOX.000000000000349935

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28 Piubelli M.L.M., Ferrulino-Schmidt M.K., Miranda R.N. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery*. JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.084

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33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery*. 2017;140(4):645-654. doi:10.1097/PRS.00000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery*. Global Open. 2021;9(3). doi:10.1097/GOX.0000000000000349935353 Establishment Labs. Dokument wewnętrzny. RTL-001020 Rationale for Motiva Implants Lifespan.

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of breast cancer at diagnosis among women with cosmetic breast implants. *Br J Cancer*. 2003; Mar 24;88(6):832-8. doi: 10.1038/sj.bjc.6600819. PMID: 12644818; PMCID: PMC2377078.

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10.21037/igs.2016.11.03. PMID: 28497021; PMCID: PMC5409903.

28. Piubelli MLM, Ferrifino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz197

29. E29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Acesso em 28 de fevereiro de 2022. <https://www.fda.gov/medical-devices/breast-implants/device-reports-breastimplant-associated-anaplastic-large-cell-lymphoma>

30. Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plast Reconstr Surg*. 2017;139(5):1042-1050. doi: 10.1097/PRS.0000000000003282

31. Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in

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a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. Journal of plastic, reconstructive & aesthetic surgery : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. Annals of surgery. 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

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34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. Plastic and Reconstructive Surgery Global Open. 2021;9(3). doi:10.1097/GOX.00000000000003499

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28 Piubelli ML, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL. ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

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25 Balk EM, Farley A, Avendano EA, Raman G. Long-Term Health Outcomes in Women With Silicone Gel Breast Implants: A Systematic Review. *Ann Intern Med.* 2 Feb 2016; 164(3):164-75. doi: 10.7326/M15-1169. Epub 10 nov 2015. PMID: 26550776.

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31 Cordeiro PG, Ghione P, NIA, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. Journal of plastic, reconstructive & aesthetic surgery : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. Annals of surgery. 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

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34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. Plastic and Reconstructive Surgery Global Open. 2021;9(3). doi:10.1097/GOX.0000000000003499.

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6 Franchelli S, Pesce M, Baldelli I, Marchese A, Santì P, De Maria A. Analysis of clinical management of infected breast implants and of factors associated to successful breast pocket salvage in infections occurring after breast reconstruction. Int J Infect Dis. 2019 Jun;71:67-72. doi: 10.1016/j.ijid.2018.03.019. Epub 2018 Apr 13. Идентификатор PubMed (PMID): 29660396.

7 Sood A, Kolamarti VS, Therattil PJ, Lee ES. Sclerotherapy for the Management of Seromas: A Systematic Review. Eplasty. 2017 Aug 28;17:e25. Идентификатор PubMed (PMID): 28890747; PMCID: PMC5575675.

8 Gioacchini M, Bottoni M, Grassetti L, Scalise A, Di Benedetto G. A simple, reliable, and inexpensive method for seroma drainage. Arch Plast Surg. 2015 May;42(3):361-2. doi: 10.5999/aps.2015.42.3.361. Epub 2015 May 14. Идентификатор PubMed (PMID): 26015895; PMCID: PMC4439599.

9 Sforza M, Husein R, Atkinson C, Zaccachedu R. Unraveling Factors Influencing Early Seroma Formation in Breast Augmentation Surgery. Aesthet Surg J. 2017 Mar 1;37(3):301-307. doi: 10.1093/asj/sjw196. Идентификатор PubMed (PMID): 28207027.

10 FDA. Risks and Complications of Breast Implants. Published 2020.

<https://www.fda.gov/medical-devices/breast-implants/risks-and-complications-breast-implants>

11 Hillard C, Fowler JD, Barta R, Cunningham B. Silicone breast implant rupture: a review. Gland Surg. 2017 Apr;6(2):163-168. doi: 10.21037/gts.2016.09.12. Идентификатор PubMed (PMID): 28497020; PMCID: PMC5409893.

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26. Hölmich LR, Lipworth L, McLaughlin JK, Friis S. Breast implant rupture and connective tissue disease: A review of the literature. *Plast Reconstr Surg*. 2007;120(7 SUPPL):1162-69. doi: 10.1097/PRS.0b00000000000003282

27. Clemens MW, Nava MB, Rocco N, Miranda RN. Understanding rare adverse sequelae of breast implants: anaplastic large cell lymphoma, late seromas, and double capsules. *Gland Surg*. 2017 Apr;6(2):169-184. doi: 10.21037/gls.2016.11.03. Идентификатор PubMed (PMID): 28497021. PMID: PMC5409903.

28. Phubelli M.M, Ferrulino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) Is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz1972 E

29. Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Докр. опубликован 28 февраля 2022 года. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breast-implant-associated-anaplastic-large-cell-lymphoma>

30. Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plastic and Reconstructive Surgery*. 2017;139(5):1042-1050. doi:10.1097/PRS.00000000000003282

31. Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery*. JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32. Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence. *Determining an Accurate Risk*. *Annals of surgery*. 2020;272(3):403-409. doi:10.1097/SLA.00000000000004179

33. Loch-Wilkinson A, Beath KJ, Knight R, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and Reconstructive Surgery*. 2017;140(4):645-654. doi:10.1097/PRS.00000000000003654

34. Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi:10.1097/GOX.00000000000003499

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²⁸ Plubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) Is Distinct from Systemic ALCL/ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. Dátum vydania: 10.1093/asj/sjz19729 E

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³¹ Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery.* JPRAS. 2020;73(5):841-846. Dátum vydania:10.1016/j.bjps.2019.11.064

³² Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. Dátum vydania:10.1097/SLA.0000000000004179

³³ Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery.* 2017;140(4):645-654. Dátum vydania: 10.1097/PRS.0000000000003654

³⁴ Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open.* 2021;9(3). Dátum vydania:10.1097/GOX.0000000000003499

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery* : JPRAS. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery*. 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery*. 2017;140(4):645-654. doi:10.1097/PRS.0000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Breast Subgroups Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi:10.1097/GOX.00000000000003499

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²⁸ Piubelli MLM, Ferrufino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sj191729 E

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³¹ Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery: JPRAS*. 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

³² Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery*. 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

³³ Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and reconstructive surgery*. 2017;140(4):645-654. doi:10.1097/PRS.0000000000003854

³⁴ Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi:10.1097/GOX.00000000000003499

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28 Piubelli MLM, Ferrulino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) is Distinct from Systemic ALCL ALCL Negative in a Patient with Gluteal Implants. *Aesthetic Surg J.* 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

29 Medical Device Reports of Breast Implant-Associated Anaplastic Large Cell Lymphoma | FDA. Accessed February 28, 2022. <https://www.fda.gov/medical-devices/breast-implants/medical-device-reports-breast-implant-associated-anaplastic-large-cell-lymphoma>

30 Doren EL, Miranda RN, Selber JC, et al. U.S. Epidemiology of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Plastic and reconstructive surgery.* 2017;139(5):1042-1050. doi: 10.1097/PRS.0000000000003282

31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after reconstruction with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery : JPRAS.* 2020;73(5):841-846. doi:10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery.* 2020;272(3):403-409. doi:10.1097/SLA.0000000000004179

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34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open.* 2021;9(3). doi:10.1097/GOX.0000000000003499

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14 Breast Cancer. Breast Implant Extrusion. Published 2020. [https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after, This is called "extrusion.](https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after, This is called)"

15 Jong, Justin MD.; Gabriel, Allen MD, FACS†; Treklel, Melissa MD.; Lawser, Amy S. MSN, RN§; Heidei, Eric PhD; Buchanan, Dallas MD, FACS; Chun, Joseph T. MD, FACS** Cohesive Round Implants and the Risk of Implant, Plastic and Reconstructive Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321 doi: 10.1097/GOX.0000000000003321

16 Khan Chopra, MD, Arvind U. Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, W. Grant Stevens, MD, Techniques to Repair Implant Malposition after Breast Augmentation: A Review, Aesthetic Surgery Journal, Volume 36, Issue 6, June 2016, Pages 660–671, <https://doi.org/10.1093/asj/sjv261>

17 Breast Asymmetry (2019). Accessed on March 26, 2021.

<https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/asymmetry>

18 Frame J. The waterfall effect in breast augmentation. Gland Surg. 2017 Apr;8(2):193-202. doi: 10.21037/gs.2016.10.01. PMID: 28497023; PMCID: PMC5409900

19 Lee Y, Song SE, Yoon ES, Bae JW, Jung SP. Extensive silicone lymphadenopathy after breast implant insertion mimicking malignant lymphadenopathy. Ann Surg Treat Res. 2017 Dec;93(6):331-335. doi: 10.4174/ast.2017.93.6.331. Epub 2017 Dec 1. PMID: 29250513; PMCID: PMC5729128.

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14 Breast Cancer. Breast Implant Extrusion. Published 2020. [https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after,This is called "extrusion."](https://www.breastcancer.org/treatment/surgery/reconstruction/corrective/implant-extrusion#:~:text=In relatively rare cases after,This is called)

15 Jong, Justin MD.; Gabriel, Allen MD.; FACS†; Trekkell, Melissa MD.; Lawser, Amy S. MSN, RNS; Heidel, Eric PhD; Buchanan, Dallas MD.; FACS; Chun, Joseph T. MD. FACS** Cohesive Round Implants and the Risk of Implant Plastic and Reconstruction Surgery - Global Open: December 2020 - Volume 8 - Issue 12 - p e3321; doi: 10.1097/GOX.00000000000003321

16 Karan Chopra, MD, Arvind U. Gowda, MD, Edwin Kwon, MD, Michelle Eagan, MD, Grant Stevens, MD, Techniques to Repair Implant Malposition after Breast Augmentation: A Review. Aesthetic Surgery Journal, Volume 36, Issue 6, June 2016, Pages 660–671. <https://doi.org/10.1093/asj/sjv261>

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28 Pibullli MLM, Ferrulino-Schmidt MC, Miranda RN. Gluteal Implant-Associated Anaplastic Large Cell Lymphoma (ALCL) Is Distinct from Systemic ALCL ALK Negative in a Patient with Gluteal Implants. *Aesthetic Surg J*. 2019;39(10):NP441-NP442. doi: 10.1093/asj/sjz19729 E

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31 Cordeiro PG, Ghione P, Ni A, et al. Risk of breast implant associated anaplastic large cell lymphoma (BIA-ALCL) in a cohort of 3546 women prospectively followed long term after treatment with textured breast implants. *Journal of plastic, reconstructive & aesthetic surgery: JPRAS*. 2020;73(5):841-846. doi: 10.1016/j.bjps.2019.11.064

32 Nelson JA, Dabic S, Mehrara BJ, et al. Breast Implant-associated Anaplastic Large Cell Lymphoma Incidence: Determining an Accurate Risk. *Annals of surgery*. 2020;272(3):403-409. doi: 10.1097/SLA.00000000000004179

33 Loch-Wilkinson A, Beath KJ, Knight RJW, et al. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plastic and Reconstructive Surgery*. 2017;140(4):645-654. doi: 10.1097/PRS.00000000000003654

34 Matros E, Shamsunder MG, Rubenstein RN, et al. Textured and Smooth Implant Use Reported in the Tracking Operations and Outcomes for Plastic Surgeons Database: Epidemiologic Implications for BIA-ALCL. *Plastic and Reconstructive Surgery Global Open*. 2021;9(3). doi: 10.1097/GOX.00000000000003499

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- | | |
|-------------------------------|----------------------------------|
| 1 = very satisfied | 4 = somewhat dissatisfied |
| 2 = somewhat satisfied | 5 = very dissatisfied |
| 3 = neutral | Leave blank if not applicable |

PACKAGING	1	2	3	4	5
1. How satisfied are you with the presentation of the external product packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you find the packaging easy to open?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the presentation of the internal packaging (thermoforms)?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Do you find the device easy to retrieve from its packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. Please suggest any packaging improvements:					
.....					
PATIENT ID CARD	1	2	3	4	5
1. How satisfied are you with the presentation of the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you think the Patient ID Card is easy to find inside the box?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the information collected on the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Please suggest any Patient ID Card improvements:					
.....					
PRODUCT INSERT	1	2	3	4	5
1. How satisfied are you with the presentation of the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfied are you with the information provided in the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Insert improvements:					
.....					

LABELING	1	2	3	4	5
1. Rate your satisfaction with the information provided on the external packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Rate your satisfaction with the information provided on the internal packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Labeling improvements:					
.....					
PRODUCT	1	2	3	4	5
1. Rate your satisfaction with the product range available inside Motiva Implant Matrix®.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfying is the appearance of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfying is the feel of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Rate your satisfaction with the implant ease of insertion into the surgical pocket.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. How satisfying is the quality of the product?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6. Rate your overall satisfaction with the product.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7. Please rate Motiva Implants® to similar products on the market. <input type="radio"/> Much better <input type="radio"/> Somewhat better <input type="radio"/> Average <input type="radio"/> Somewhat worse <input type="radio"/> Much worse If your answer is "somewhat worse", or "much worse", please explain:					
8. Would you use Motiva Implants® again in the future? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain:					
9. Would you recommend Motiva Implants® to a colleague? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain:					
SUPPORT	1	2	3	4	5
How satisfied are you with delivery of the products?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
How satisfied are you with the support provided by your distributor?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Please suggest how Establishment Labs can improve your customer experience:					
.....					

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MOTIVA IMPLANTS FEEDBACK FORM

This survey provides valuable customer input that helps us offer the best possible service. **Thank you** for completing and returning the Motiva Implants® Feedback Form by mail to Establishment Labs Coyol Free Zone & Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica, or by email: support@doctors-motivaimplants.zendesk.com.

Name:	Implant Information
Gender:	Reference Number:
Address:	Serial Number:
City:	Product:
Country:.....	Lot Number:

Please rate your satisfaction level with each of the following statements:

1 = very satisfied 2 = somewhat satisfied 3 = neutral	4 = somewhat dissatisfied 5 = very dissatisfied Leave blank if not applicable
--	---

PACKAGING	1	2	3	4	5
1. How satisfied are you with the presentation of the external product packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you find the packaging easy to open?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the presentation of the internal packaging (thermoforms)?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Do you find the device easy to retrieve from its packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. Please suggest any packaging improvements:					

PATIENT ID CARD	1	2	3	4	5
1. How satisfied are you with the presentation of the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you think the Patient ID Card is easy to find inside the box?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the information collected on the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Please suggest any Patient ID Card improvements:					

PRODUCT INSERT	1	2	3	4	5
1. How satisfied are you with the presentation of the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfied are you with the information provided in the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Insert improvements:.....					

LABELING	1	2	3	4	5
1. Rate your satisfaction with the information provided on the external packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Rate your satisfaction with the information provided on the internal packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Labeling improvements:					
PRODUCT	1	2	3	4	5
1. Rate your satisfaction with the product range available inside Motiva Implants®.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfying is the appearance of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfying is the feel of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Rate your satisfaction with the implant ease of insertion into the surgical pocket.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. How satisfying is the quality of the product?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6. Rate your overall satisfaction with the product.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7. Please rate Motiva Implants® to similar products on the market. <input type="radio"/> Much better <input type="radio"/> Somewhat better <input type="radio"/> Average <input type="radio"/> Somewhat worse <input type="radio"/> Much worse If your answer is "somewhat worse", or "much worse", please explain:					
8. Would you use Motiva Implants® again in the future? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain:					
9. Would you recommend Motiva breast Implants® to a colleague? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain:					
SUPPORT	1	2	3	4	5
How satisfied are you with delivery of the products?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
How satisfied are you with the support provided by your distributor?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Please suggest how Establishment Labs can improve your customer experience:					

MOTIVA IMPLANTS® FEEDBACK FORM

This survey provides valuable customer input that helps us offer the best possible service. Thank you for completing and returning the Motiva Implants® Feedback Form by mail to Establishment Labs Coyal Free Zone & Business Park Building, 4th Street, Building B-15, Alajuela, Costa Rica, or by email: support@doctors-motivaimplants.zendesk.com

Name:.....	Implant Information
Gender:..... Age:.....	Reference Number:.....
Address:.....	Serial Number:.....
City:.....	Product:.....
Country:.....	Lot Number:.....

Please rate your satisfaction level with each of the following statements:

- | | |
|------------------------|-------------------------------|
| 1 = very satisfied | 4 = somewhat dissatisfied |
| 2 = somewhat satisfied | 5 = very dissatisfied |
| 3 = neutral | Leave blank if not applicable |

PACKAGING	1	2	3	4	5
1. How satisfied are you with the presentation of the external product packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you find the packaging easy to open?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the presentation of the internal packaging (thermoforms)?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Do you find the device easy to retrieve from its packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. Please suggest any packaging improvements:					

PATIENT ID CARD	1	2	3	4	5
1. How satisfied are you with the presentation of the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you think the Patient ID Card is easy to find inside the box?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the information collected on the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Please suggest any Patient ID Card improvements:					

PRODUCT INSERT	1	2	3	4	5
1. How satisfied are you with the presentation of the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfied are you with the information provided in the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Insert improvements:					

LABELING	1	2	3	4	5
1. Rate your satisfaction with the information provided on the external packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Rate your satisfaction with the information provided on the internal packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Labeling improvements:					
PRODUCT	1	2	3	4	5
1. Rate your satisfaction with the product range available inside Motiva Implant Matrix®.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfying is the appearance of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfying is the feel of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Rate your satisfaction with the implant ease of insertion into the surgical pocket.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. How satisfying is the quality of the product?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6. Rate your overall satisfaction with the product.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7. Please rate Motiva Implants® to similar products on the market. <input type="radio"/> Much better <input type="radio"/> Somewhat better <input type="radio"/> Average <input type="radio"/> Somewhat worse <input type="radio"/> Much worse If your answer is "somewhat worse", or "much worse", please explain:					
8. Would you use Motiva Implants® again in the future? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain:					
9. Would you recommend Motiva Implants® to a colleague? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain:					
SUPPORT	1	2	3	4	5
How satisfied are you with delivery of the products?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
How satisfied are you with the support provided by your distributor?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Please suggest how Establishment Labs can improve your customer experience:					

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5. Please suggest any packaging improvements:

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Copyright ©2023 Establishment Labs S.A.